About The User Manual

Thank you for purchasing our product!

In order to enable you to skillfully operate our product as soon as possible, a detailed user manual is attached.

Please make sure to read all the content when installing and using the product for the first time.


To improve the performance and reliability of its parts, the product (including hardware and software) may be changed from time to time, during which, we will try to modify or add contents. Please forgive us as there may still be inconsistency with some descriptions.

We look forward to your corrections in case of any errors and omissions in this manual.



All rights reserved © Shenzhen Biocare Bio-Medical Equipment Co., Ltd.

Information contained in this document is proprietary to Shenzhen Biocare Bio-Medical Equipment Co., Ltd. No part of this publication may be reproduced, stored in a retrieval system, or transmitted in any form, in whole or in part, by any means electronic, mechanical or other, including photocopying and recording, for any purpose without written permission of Biocare.





The CE mark is a protected conformity mark of European

Community. The products herewith comply with the requirements of the

Medical Device Directive 93/42/EEC



Shanghai International Holding Corp. GmbH (Europe)

Eiffestraße 80

20537 Hamburg

GERMANY



Shenzhen Biocare Bio-Medical Equipment Co., Ltd.

#A735,7/F, Block A, Shenzhen Mingyou Industrial Products

Exhibition&Procurement Center, Baoyuan Road, Xixiang Sub-district,Bao'an

District, 518102 Shenzhen, PEOPLE'S REPUBLIC OF CHINA

Tel: 86 -755 -36615333 Fax: 86-755-27960643

Website: http://www.biocare.com.cn

CAUTION: In US, Federal law restricts this device to sale by or on the order of a physician. Please read the user manual carefully prior to use.

Explanation of key words

•* WARNING
Indicate a potentially hazardous situation which, if not avoided, could result in death or serious injury.
A CAUTION
Indicate a potentially hazardous situation which, if not avoided, may result in slight personal injury or equipme
failure.
EXPLANATION
Other important information besides warning or caution.

Explanation of part symbols

Symbol	Explanation	Symbol	Explanation
~	Alternating Current		Type CF Applied Part
CB	Direct Current	- ** -	Type CF applied part including defibrillation-Proof
→ □	Battery Charging	몲	LAN Port
4	Equipotential Terminal	←	USB Port
\triangle	Please refer to user manual!		

Conventions

Format	Explanation	
******	Used to quote the texts in the screen of the machine.	
[****]	Used to quote the shortcut buttons or keys in the screen of the machine.	
TEXT Used to quote the referenced chapters or sections in this manual.		

Content

Content	I
Foreword	1
Chapter 1 Introduction	9
1.1 Equipment Overview	10
1.1.1 Front View	10
1.1.2 Left View	11
1.1.3 Right View	11
1.1.4 Operating Buttons	12
1.2 Waveform Display	13
1.2.1 Same Screen Display	13
1.2.2 Split-screen Display	14
Chapter 2 Preparation	15
2.1 Locate the ECG Machine	15
2.2 Install Recording Paper	15
2.3 Connect to Power Supply	16
2.4 Connect to Patient Cable	17
2.5 Power On/Off	17
2.6 Connect to Network	18
2.7 Apply Electrodes	19
2.7.1 Electrodes Attachment	20
2.7.2 Lead Signals Formation Scheme	21
Chapter 3 Entering Patient Information	23
3.1 Enter Patient Information	23
3.2 Introduction of Input Method	25
3.2.1 Standard Character Keyboard	25
3.2.2 Digital Keyboard	26
Chapter 4 ECG Recording	27
4.1 Introduction of Sensitivities, Filters, Print Speed	28

4.2 Record ECGs	29
4.2.1 Main Steps to Record ECG	30
4.3 Introduction of Record Mode	31
4.3.1 Auto Mode	32
4.3.2 Manual Mode	
4.3.3 Upload Mode	34
4.4 Advanced Mode	34
4.4.1 Rhythm Mode	34
4.4.2 Analysis Report Mode	35
Chapter 5 Setting System Parameters	37
5.1 ECG Setting	37
5.2 Print Setting	39
5.3 Display Setting	40
5.4 Patient Information Setting	40
5.5 System Setting	40
5.6 Factory Maintain	42
Chapter 6 Data Management	43
6.1 Open an ECG File	44
6.2 Edit an ECG File	44
6.3 Delete ECG Files	44
6.4 Copy and Move ECG files	45
Chapter 7 Maintenance	
7.1 Main Unit	
7.2 Patient Cable	47
7.3 Cleaning and Disinfection	
7.4 Recording Paper	
7.5 Battery	
7.6 Silicon Rubber Shaft for Printing	
7.7 Thermal Print Head	
Chapter 8 Troubleshooting	
- ··· t · · ·	

8.1 Lead Fault	51
8.2 Printer Failure	52
8.3 Indicator of Lead Off	52
8.4 AC Interference	53
8.5 EMG Interference	53
8.6 Baseline Wander	54
8.7 The ECG Machine cannot be turned on	54
8.8 Paper Feeding Failure	55
8.9 Battery is quickly Charged and Discharged	55
8.10 Wrong Analysis Result	55
8.11 File Uploading Failure	56
Appendix A Package and Accessories	57
A.1 Packing List	57
A.2 Dimensions and Weight	57
Appendix B Technical Specification	59
B.1 Specifications	59
B.1.1 Main Unit	59
B.1.2 Recorder Specification	60
B.1.3 Other Specification.	60
B.2 Environment Requirements	61
Appendix C Working Principle and Block Diagram	63
C.1 Power Supply Subsystem	63
C.2 Acquisition Module	64
C.3 Control System	65
Appendix D List of Interpretation Codes and Corresponding Description	67
Appendix E Measurement, Diagnosis, Analysis and Assessment of ECG Machine	75
E.1 Methods to determine the amplitude of P, QRS, ST and T wave	75
E.2 Processing method of isoelectric segment in QRS complex	76
E.3 Low incidence heart disease not included in testing and diagnosis database	76
E.4 ECG diagnosis categories and the number of ECG test of each category	77

	E.5 The smallest waveform identified by the device and the stability of measurement when noise exists 77
	E.6 Low incidence cardiac rhythm not included in the ECG rhythm test database
	E.7 ECG rhythm diagnosis categories and ECG test number of each category
	E.8 Sensitivity regularly test instructions
	E.9 Distortion test
App	pendix F EMC-Guidance and manufacture's declaration81
	F.1 Guidance and manufacturer's declaration – electromagnetic emission
	F.2 Guidance and manufacture's declaration-electromagnetic immunity for all EQUIPMENT and SYSTEMS
	F.3 Guidance and manufacture's declaration-electromagnetic immunity for EQUIPMENT and SYSTEMS
	those are not LIFE-SUPPORTING
	F.4 Recommended separation distance between portable and mobile RF communications equipment and the
	EQUIPMENT or SYSTEM for EQUIPMENT or SYSTEM that are not LIFE-SUPPORTING

Foreword

A CAUTION

- This ECG machine shall be used by qualified health professionals in the medical units, who need to analyze the ECG waveforms and give diagnostic results.
- In sequence to use this ECG machine correctly, safely and effectively, please read through the user manual carefully.

★ Safety Information

WARNING

- Avoid using and storing in the places with sulfur, salt, alkaline gas or with risk of gas leakage.
- Avoid using in the places with anesthetic gases, flammable gases such as oxygen, hydrogen or other flammable chemicals, or it may cause explosion or fire.
- Select a room with intact infrastructure (good power supply system and grounding facilities).
- Be cautious when the patient is connected with more than one instrument, because the total leak current may be harmful to the patient. Devices in compliance with the standard of IEC60601-1 are allowed to be connected to this ECG machine, and the equipotential points of all the connected devices should be connected reliability. (The equipotential point and the protection ground of this ECG machine have been connected). Total leak current should be measured by the users to determine that it meets the requirement and can be used after connection.
- All the analog and digital equipment that is connected to this ECG machine in the patient environment has to be in compliance with the standard IEC60601-1; All the analog and digital equipment that is connected to this ECG machine out of patient environment has to be in compliance with other national safety standards (IEC or ISO safety standards); the composition system should be in compliance with the standards of IEC 60601-1-1.
- When the equipment is used simultaneously with cardiac defibrillators, avoid contacting with patients or hospital beds. All the electrodes connected and unconnected to patients as well as patients themselves do not have to be grounded. Do not use other electrical stimulators at the same time. If needed, it should be a professional technician to carry out the operation.

- Chest and limb electrodes along with the device in the packing box could not meet the requirements of defibrillation polarization recovery time (however, they are essential accessories of ECG), should not be used immediately for reliable measurements and diagnostics after defibrillation. To ensure proper defibrillator protection, use only the recommended disposable electrodes (Name: Skintact; Type, RT-34), lead wire and electrode adapters by our company. To ensure the protection of defibrillator discharge, use the lead wire with defibrillation-proof by our company.
- When the ECG machine is used together with a defibrillator or other electrical stimulators (like high-frequency surgical devices), we recommend using disposable chest electrodes. Otherwise the patient may get serious injury by using mental electrodes.
- During defibrillation, the device can detect the discharge of defibrillator, and process automatically, and then quickly recover the waveform.
- Keep the ECG machine electrodes away from the electrodes of high-frequency electrosurgical units. Ensure the resistance between the electrosurgical unit and patient body is as low as possible. If necessary, the disposable electrodes can be used because of its larger contact area on the human body, which can keep the high-frequency current density in an acceptable range.
- When the relevant packaging material including, depleted batteries and scrapping products are disposed, please follow the local laws; the user should properly follow requirements of local laws, and recycling laws.

A CAUTION

- Avoid contacting with water or other liquids, and avoid using and storing in spaces with too large barometric pressure, humidity and temperature beyond the prescribed standards, poor ventilation, or with excessive dust.
- The ECG machine should be place on the flat horizontal table and avoid strong vibration and mechanical shock while moving.
- The frequency of AC power supply and system voltage should be complied with the requirements. More importantly, the current capacity should be sufficient.
- The instrument should not be surrounded by high-voltage cables, ultrasound equipment, electrotherapy machines and other high-power equipment.
- To more accurately record the ECG, the equipment shall be placed in a quiet and comfortable environment.
- The circuit of application parts works based on floating ground and meets the safety standards in IEC60601-1 CF Type. It can be used in acquiring the body surface ECG signals, but cannot applied to the heart directly.

- Turn off the ECG machine if an accident happens.
- Please clean and disinfect reusable electrodes with medical alcohol before using.
- The conductive parts of electrodes and connectors (including neutral electrodes) on the ECG machine should not in touch with other conductive parts.
- Do not press the buttons with sharp or hard articles or it may cause permanent damage to the buttons.
- Do not make any modifications to this ECG machine.
- Perform regular maintenance and inspection for this ECG machine and all its accessories (at least once every six months).
- The maintenance and repair of this ECG machine should be performed by experienced technicians. When there is any functional abnormality, it should be clearly identified to prevent the ECG machine from running with fault.
- The electrical schematic diagrams and parts listed are only provided to a qualified repair station or technicians recognized by the company.

★ General Operation Precautions

Before operation:

≜ CAUTION

- Make sure the ECG machine is in good condition and the recording paper is sufficient.
- Make sure the temperature and humidity of operating environment comply with the requirements.
- Do not operate ECG machine in a place with X-ray equipment, ultrasound scanner, or other similar equipment. Those equipments may interfere with the ECG machine. If necessary, power off the mentioned equipments or move the ECG machine to an environment without interference.
- Make sure all lead wires are connected correctly and are kept away from the AC power cable.
- Make sure the equipotential cable of the ECG machine is reliable and properly connected.
- Make sure the power cable is properly connected with the ECG machine and is not twisted with other cables or wires.
- Put the lead wires in good order before connecting them to the electrodes.
- Make sure the electrodes are in good contact with skin. Please refer to *Apply Electrodes* for details.
- Please install the ECG machine near an AC power outlet. Cut off the power supply immediately when there is an emergency.

- If the patient is nervous, please explain to the patient that ECG examination is easy and no harm.
- Please keep the patient quiet and motionless.
- Use wide hospital beds and keep the patient from touching the metal parts of the hospital bed, which may cause interference in ECG waveform recording.
- Keep the Examination room quiet and comfortable.

WARNING

- All circuits that come in contacting the patient directly should be examined closely.
- When using the battery as the power supply, please check the voltage and condition of the battery first and make sure the battery is fully charged. For new batteries, please discharge and charge it fully before using.
- Use only 3-core power cable when using AC power, otherwise hazard of electric shock to the patient and operator cannot be completely eliminated. If the power cable is not working, only the built-in battery can safely power the ECG machine.
- Make sure equipotential connection is complete and reliable, or else only use the built-in battery.

In operation:

WARNING

- The physician should observe the patients closely without leaving during the operation. If necessary, turn off the ECG machine and remove the electrodes to ensure patient's safety.
- Prevent the patients from contacting the other parts of the ECG machine or other conductors except for the electrodes.

After operation:

≜ CAUTION

- Please return to the main interface before turning off the ECG machine.
- Remove the electrodes gently and do not pull the lead wires emphatically.
- Clear up the ECG machine and all the accessories for trouble-free operation of next use.

About LCD screen

CAUTION

- Do not place any heavy objects onto the LCD screen or strike it, otherwise it could damage the LCD screen.
- When not using it, please put it away or have a cover on it. Keep it away from water.

About lithium battery

WARNING

- Only the authorized installation or service engineer can open the battery cover and replace the battery; do use the same type of rechargeable lithium battery provided by our company.
- The positive and negative terminals of the batteries cannot be reversed, or it could cause an explosion.
- Do not connect the two polarities of the battery with metal wires. Otherwise, there will be the hazard of fire.
- Do not use the battery near a heat source or in and the environment with temperature up to 60 °C; do not heat the battery or throw it into the fire.
- Keep the battery away from water; do not drop the battery into the water.
- Do not press any metal into the battery; Do not hammer or beat the battery or use other ways to damage the battery, otherwise it will cause heat, smoke, deform or burn, which is very dangerous.
- When you find battery leakage or its emitting unpleasant odors, please get away from it immediately. If the fluid leaks onto the skin or clothes, wash with clean water at once. If the electrolyte enters the eyes, do not rub the eyes, wash with clean water immediately, and then go to see a doctor.
- The users need to check the battery status regularly. When the battery reaches the end of its lifetime, when it smells, deforms, discolor, contorts, the users should stop using and dispose of it according to local regulations.

★ EMC Considerations

This ECG machine conforms to the IEC60601-1-2, a safety standard for medical electronic devices or systems. However, the electromagnetic environment exceeding the limit or level defined by the standard IEC60601-1-2 will introduce the unwanted interference to the ECG machine, disable its intended functions or it will compromise its intended performance. Thus, if there is any discrepancy with this ECG machine compared to its intended

functions during operation, please do not use it for longer until the adverse affect is identified and eliminated. The appropriate preventing measures are given below by this manual for such cases:

■ Influence of radiated electromagnetic wave:

The use of a mobile phone may affect this ECG machine. Instruct all the people around to turn off their mobile phone or mini-radio devices when any medical electronic device is in use.

■ Influence of impact and conductive electromagnetic waves:

The high frequency noise produced by other devices can be introduced into this ECG machine through the alternating current socket. Please identify the noise source first, and if possible, stop the working of related devices. If they are not allowed to be stopped, measures such as application of noise abatement device should be taken to minimize the influence.

■ Influence of static electricity:

The static electricity in a dry environment (indoor) may affect this ECG machine, especially in winter. Please humidify the indoor air or pre-discharge the static electricity on the cable and the electrocardiogram recording personnel prior to using this ECG machine.

■ Influence of thunder and lightning:

A thunder and lightning strike nearby may cause voltage surge in this ECG machine. You can unplug the power supply and run the ECG machine using its internal battery in case of any danger.

Instrument classification

Methods	Class
By Type of Protection Against Electric Shock	Class I, internal power supply
By Degree of Protection Against Electric Shock	Type CF applied part
By Degree of Liquid Proof	Ordinary equipment (enclosed device without liquid proof)
	This equipment is unsuitable for use in an environment with
By Level of Protection Against Explosion	air, oxygen or nitrous oxide mixed with flammable anesthetic
	gas.
By Mode of Operation	Continuous operation equipment

Maintenance Warranty

Our company guarantees the new instrument on the material and technological qualification for this product within 18 months and the accessories within 6 months since purchasing day, while consumables are not covered by the guarantee in principle. This guarantee is also inapplicable to the products undergoing any modification, disassembly, refitting or self-repairing without permission of our company, as well as the products damaged by accidents, fire disaster, thunder and lightning, flood and other disaster, intentional damage, improper installation and improper usage.

A CAUTION

- For all dated reference documents in this manual, its subsequent amendments (excluding corrections) or revisions do not apply to this manual; for undated reference documents, the latest version applies in this manual.
- Due to product improvement, the content of this User Manual may differ from the product you purchased, which will not affect the usage, please operate according to the actual functions of the product.
- This manual is subject to change without prior notice. We apologize for any inconvenience caused.

-- Blank Page --

Chapter 1 Introduction

[Main structure, performance]

The ECG machine is mainly composed of the main unit, lead wire, limb electrodes and chest electrodes.

[Scope of application]

This ECG machine is used to extract the electrocardio complex from the human body for clinical diagnosis and research.

[Intended Use]

- The diagnostic applications include: check the cardiac abnormalities of the general population; detect the chest pain in patients with acute myocardial ischemia and myocardial infarction, and evaluate the patients with arrhythmias;
- Suitable for: adults (older than 12 years old), pediatrics (age between 29 days to 12 years old), and neonates (infants born less than 28 days after 37 weeks to 44 weeks of pregnancy);
- Used in: hospitals, clinics;
- The automatic analysis program of this ECG machine focuses on the high sensitivity of detecting high-risk patients with cardiac abnormalities.

1.1 Equipment Overview

1.1.1 Front View




Figure 1-1 Front View

Number	Name	Description
1	Switch Button	Push down to open the Paper Drawer.
2	Paper Drawer	Place the recording paper.
3	Display Screen	Display the waveforms, patients' information and the device status.
	Operating Buttons	For button operations and inputting methods. Refer to <i>Operating</i>
4		Buttons for more details.

1.1.2 Left View




Figure 1-2 Back View

Number	Name	Description
1	Power Supply Socket	Connect to the AC power adaptor.
2	Switch Button	Push down to open the Paper Drawer.
2	Equipotential Terminal	Connect to the equipotential cable.

1.1.3 Right View




Figure 1-3 Side View

Number	Name	Description
1	ECG Patient Cable Connector	Connect to the Patient Cable.
2	USB Ports	Insert the USB device to save data; Insert the Bar Scanner.
3	LAN Port	Connect to net cable.

1.1.4 Operating Buttons




Figure 1-4 Keyboard

Number	Name	Description
1	6/6)	Press to power on or off the ECG machine.
2	Ŷ?	Press to set patient information.
3		Indicate the power status of the ECG machine. From left to right: AC power indicator, battery indicator and
3	~ / 🝱 / 🗀	charging indicator.
4	(1mV/GP)	In Manual Mode, press to print the calibration waveforms of 1mV to check the current sensitivity.
5	(((((((((((((In Auto Mode, press to copy the previous report. Press to start or stop printing the ECG waveforms and report.
6&7	(ENTER)	Press Up/Down/Right/Left button to select a menu or an option. Press [ENTER] to confirm, open a submenu, or toggle between two options in a submenu.
8		In Main interface, press to enter configuration Menu. In other interface, press to exit.
9		Press to select a sensitivity.

10	2	Press to set low-pass filter, baseline wander filter and AC filter.
11		Press to select a record mode and record format.
12	(Lead) (Lead>)	In manual mode, press to switch among different lead groups. In split-screen mode, press to switch among different screens.

1.2 Waveform Display

EXPLANATION

Screen display may slightly differ from the product you purchased, which will not affect your usage. Please operate according to the actual functions of the product.

In same screen displayed interface, 12 lead waveforms will be displayed on one interface.

In split-screen displayed interface, 12 lead waveforms will be displayed on several interfaces, which make it possible to show the waveform details more clearly.

Select [Display], set the display format and lead format.

1.2.1 Same Screen Display

Waveform display in same screen, 3 X 4 lead format:




Figure 1-5 Main Interface

Waveforms display in same screen,

6 X 2 lead format:

ī	V1	
II	V2	
	V2 V4	
aVR	V3	
aVL	V5	
aVF	V6	

Waveforms display in same screen,

12 X 1 lead format:

1			
II			
Ш			
aVR			
aVL			
aVF			
V1			
V2			
V3			
V4			
V5			
V6			

1.2.2 Split-screen Display

Waveform under split- screen, 3 X 4 lead format:




Figure 1-6 Lead Switch

Press (Lead) to change displayed waveforms.

Chapter 2 Preparation

2.1 Locate the ECG Machine

Please refer to Foreword.

2.2 Install Recording Paper

EXPLANATION

Installation of the recording paper may slightly differ from the product you purchased, which will not affect your usage, please operate according to the actual feature of the product.

See the descriptions below to install the rolling paper:




Figure 2-1 Install the Recording Paper

- 1. Push the button downwards.
- 2. Open the Paper Drawer cover.

- 3. Remove the roller.
- 4. Insert the roller into the Rolling Paper.
- Load the recording paper fitly into the Paper Drawer. And make sure the grid side of the paper facing downwards and towards the Thermal Print Head.
- 6. Pull out the paper about 2cm and press to close the Paper Drawer cover.

A CAUTION

- Please make sure the recording paper is installed fitly and straightly, otherwise it may be stuck.
- If paper is absent, or used up, or not placed well, alarm will appear on the main interface and the machine does not print.

2.3 Connect to Power Supply




Figure 2- 2 Connect to AC Power Supply

- 1. Connect the AC Power Cable to the ECG machine and the Power Outlet
- 2. Connect Equipotential Cable to the ECG machine and the Equipotential Terminal in the room.

EXPLANATION

The ECG machine is equipped with a built-in rechargeable battery and requires no extra installation. The operator needs to check the battery's capacity before usage.

CAUTION

■ When the ECG machine is operated together with other medical equipment, please use the accompanying -- 16 -- User Manual for Electrocardiograph

Equipotential Cable and connect the Equipotential Terminal of the ECG machine altogether with that of the other equipments so as to protect the patient from possible electric shock due to current leakage from those equipments.

■ It is forbidden to connect the Equipotential Cable to a conductive water pipe or other pipes. Otherwise, there may be hazard of electric shock to the patient.

2.4 Connect to Patient Cable




Figure 2-3 Connect to Patient Cable

Connect the Patient Cable to the ECG machine.

A CAUTION

Do not use any other Patient Cable except that supplied one. The ECG Patient Cable Connector is exclusively used for connecting the Patient Cable. Do not use it for other purposes.

2.5 Power On/Off



to power on/off the ECG machine.

The ECG machine enters Standby mode if it is not in use for a setting period. Set the duration in [System Setting] > [Standby Time]. To exit standby mode, press any key.

The ECG machine will shut down automatically if it is not in use for a setting period. Set the duration in [System Setting] > [Auto Power Off].

2.6 Connect to Network

A CAUTION

In the data transmission, if the ECG machine warns "Network connection failed", please reconnect to network or reset the network configuration.

- 1. As shown in the following figure, the cable network system is composed of the ECG machine, the switchboard and the server.
- 2. Go to [Menu] > [System Setting] > [Cable Network] and set the [IP address], [Subnet Mask] and [Default Gateway] of the ECG machine. If the IP address is in the same network segment of the server, the subnet mask and gateway shall be as the set value of the server. If the IP address is not within the same network segment of the server, subnet mask and gateway of the [Cable Network] shall be set according to actual situations. And make sure the specified gateway does support the data transmission between the two network segments.
- 3. Select [Server Setting] and set the correct [IP address] and [Port] number of the server.
- 4. When the machine is connected to network successfully, the icon will display in the main interface.




Figure 2- 4 Cable Network Connection

2.7 Apply Electrodes

Before attaching the electrodes to patient, wipe the skin where the electrodes attached by using medical alcohol, and then apply ECG gel on the skin. After that, place electrodes to the right position. If using the vacuum ball electrodes, you should apply ECG gel on the electrodes, and then pinch the suction ball to make sure the electrodes are attached to skin tightly.

A CAUTION

- Proper electrode attachment is vital for obtaining accurate ECG waveforms; therefore, please ensure good contact between the skin and electrodes.
- Do not use the new electrodes and the used ones at the same time. Replace all electrodes together when any one of them is supposed to be replaced.
- Don't use disposable electrodes more than one time.
- Confirm the disposable electrodes are within the validity period.
- Use the disposable electrodes as soon as possible after opening the package (generally within 7 days).
- Electrodes or conducting point of lead wires shall not be in contacting with any other metal part or conductor.
- Avoid the electrodes to be dragged by the lead wires.
- Make sure patient's skin contacted with the electrodes has been well pretreated.
- Clean the stain on the electrodes with medical alcohol whenever the electrodes are contaminated.
- Make sure metal electrodes of limb electrodes touch fully with skin and tightly enough.
- Make sure adjacent electrodes and ECG gel, especially chest ones, are not contacted with each other.
- If the examination involves a short period of time, if ECG gel is unavailable, please wipe the skin with medical alcohol to keep the skin clean and moist, rapidly attach the electrodes.
- It's not allowed to use saline water as substitute when ECG gel is not available. The saline water will cause corrosion on the electrodes.
- For chest and back application of pediatric, disposable electrodes are suggested.
- Electrodes shall be properly stored. When electrodes have been used for a certain period, they may become corroded and oxidized at the surface. Whenever this happens, the electrodes must be replaced.
- Do not mix electrodes of different types and manufacturers. Do not use re-useable electrodes and disposable ones together, or it will affect the recording.

■ Please use our company's or authorized electrodes to make sure qualified ECG signals.

2.7.1 Electrodes Attachment

Lin	Limb Electrodes Placement				
	IEC	АНА	Description	Figure	
	R Red	RA White	Right Arm		
	L Yellow	LA Black	Left Arm		
	N Black	RL Green	Right Leg		
	F Green	LL Red	Left Leg		
Sta	ndard 12-lead	Attachment			
	IEC	AHA	Description	Figure	
A	C1 Red	V1 Red	Fourth inter-costal space at right sternal border.		
В	C2 Yellow	V2 Yellow	Fourth inter-costal space at left sternal border.		
С	C3 Green	V3 Green	Equidistant between B and D	a e o o e	
D	C4 Brown	V4 Blue	Fifth inter-costal space at left mid-clavicles line		
Е	C5 Black	V5 Orange	Left anterior auxiliary line at the horizontal level of D		
F	C6 Purple	V6 Purple	Left mid- auxiliary line at the horizontal level of D		

2.7.2 Lead Signals Formation Scheme





Figure 2- 5 Lead Connection

-- Blank Page --

Chapter 3 Entering Patient Information

3.1 Enter Patient Information



Enter the patient information

You can enter patient ID Number, patient name, gender, and age of years. Refer to *Set Patient Information* to get more details.

Input ID number:




Figure 3- 1 Inputting Patient Information

Select the arrow button and press **ENTER**, and a menu with below three options pops up.

[Automatic Coding]: ID code is automatically generated by ECG machine when admitting a new patient and an ID code will be automatically increased each time when you press ENTER.

[Manual Coding]: you can code according to his own demand with numbers and letters.

[Barcode Scanner]: you can directly scan the bar code using the scanner to generate ID code.



Figure 3-2 Barcode Scanning

EXPLANATION

When Barcode Scanner is chosen, the onscreen keyboard will be disabled. For the usage of the Barcode Scanner, please refer to the User Manual of your Barcode Scanner.

A CAUTION

- Improper patient information may lead to misdiagnose. Please check the information for each new patient carefully.
- Please avoid confusing the patients' ID numbers. Otherwise, it may cause ECG data loss or mistake.

3.2 Introduction of Input Method

3.2.1 Standard Character Keyboard

Select a text box and press **ENTER** to open the onscreen keyboard. Use Up/Down/Right/Left button to select the letters and then press **ENTER** to confirm.

Press to exit the onscreen keyboard and return to the previous page.

Or, select **End** button on the screen and press **ENTER** to exit the onscreen keyboard.




Figure 3-3 Standard Character Keyboard

Number	Name	Description	
1	Characters Area	Input letters or punctuations.	
2	Delete	Delete previous inputted character.	
3	Cursor Moving	Move the position of cursor on the screen.	
4	End	Exit the onscreen keyboard.	
5	Page Down	/	
6	Page Up		
7	Space	Input null characters.	
8	Symbols	Switch to symbols pad to input kinds of symbols.	
9	Caps	Switch between uppercase and lowercase letters inputting.	

3.2.2 Digital Keyboard

Select a text box and press **ENTER** to open the digital keyboard. Use Right/Left button to select the numbers and then press **ENTER** to confirm. Use Up/Down/ to exit the digital keyboard.

Or, press to exit the digital keyboard and return to the previous page.




Figure 3- 4 Digital Keyboard

Number	Name Description	
1	Characters Area	Input numbers or punctuations.
2	Cursor Moving	Move the position of cursor on the interface.
3	Delete	Delete previous inputted character.
4	End	Exit the Digital Keyboard.
5	Space	Input null characters.

Chapter 4 ECG Recording

After the ECG machine has been powered on and all the leads are well connected, the following main interface will be displayed. And the ECG machine is ready for recording.



Figure 4 - 1 Main Interface

Number	Name	Description	
1	Heart Rate Icon	Display patient's heart rate.	
2	Lead Status	Display the location and state of the electrode in the human body. Select for a larger view. If the lead wires are not connected properly, such as electrode falls off, the corresponding electrodes on the view will flash to alarm.	
3	Record Setting	Display record mode and record format; Use to select other record mode and format.	
4	Filter, Sensitivity and Print	Display the current Filter, Sensitivity, and Print Speed.	

	Speed Status Area			
5	Patient Information	Display patient information;		
		Display text alarm information, including: system failure about		
6	Alarm Area	Patient Cable/Print head/Paper, lead off, AC interference, EMG		
		interference, Baseline wander, and Data overflow, etc.		
7	System Status	Indicates system status, for example, mute, recording, network,		
/		USB connecting, battery, etc.		
8	Shortcut Keys	Quick operation to set up parameters and execute functions.		
9	Waveform Display	Display real time waveform.		

4.1 Introduction of Sensitivities, Filters, Print Speed

Before printing, set the following configuration:

Press to select the desired sensitivites.

Press to select the desired Low-pass Filter, Baseline Filter and AC Filter.

On the main interface, use navigation buttons to open the Print Speed Setting Menu and select the desired print speed.




Figure 4- 2 Setup of the Sensitivities, Filters, and Print Speed

EXPLANATION

- Noise may affect the quality of the ECG signals. You can change the filter parameters to optimize the displayed or printed ECG waveforms.
- The filters settings affect the displayed and printed ECG waveforms only, not the analysis result.
- In order to reduce baseline interference, a baseline wander filter should be employed. To make sure the ST segment is not distorted, the AAMI standards recommend that cut-off frequency of the baseline wander filter is lower than 0.67Hz.

4.2 Record ECGs

EXPLANATION


- Pre-acquisition only works in Auto Mode and Economic mode:

 When the Pre-acquisition mode is enabled, the ECG machine print and upload waveforms automatically; when the Pre-acquisition mode is disabled, press to print or upload waveform.
- [Record Format] is a waveform pattern traced on the recording paper. Please see *Technical Specification* for specific recording formats.
- [Synchronous] and [Real Time] are activated only when waveforms are printed in more than one column.

 [Synchronous]: all the start point of the waveforms are the same;
 - [Real Time]: the start points of the waveforms in the same columns are the same; but the start points of different column are continuous to the end point of the waveforms in the previous column.
 - Set them in [Print Setting] > [Print Data Type].
- If you have blank paper, select [**Print Setting**] > [**Print Grid**] > [**Enable**] to print the ECG showing the grids.

 If you have paper with grids, suggest selecting [**Disable**] to print the ECG without showing the grids.


4.2.1 Main Steps to Record ECG



Set up the record mode and record format



Observe the waveforms display



Exit the recording

When the waveforms recording is about to finish (under any record mode except manual mode), the system will carry on the interpretation automatically. Please refer to *Auto Mode*.

A CAUTION

- After the heart rate and waveforms are stable, you can print out the ECG waveforms together with the interpretation.
- This ECG machine can detect the lead connecting status continuously, and if leadoff is detected, the corresponding lead code will flash in the Alarm Area on the main interface, accompanied by sound alarm. When "Lead off" continues, please check carefully the connections from skin to the ECG machine (including electrodes and lead wires). Alarm will disappear when connections become reliable.


4.3 Introduction of Record Mode

EXPLANATION

- Set up the Waveform length to be printed and uploaded in [Menu] > [ECG Setting] > [Wave Sample Time].
- In [ECG Setting], if [Auto Upload] is enabled, the ECG machine will automatically upload the waveforms and the reports after printing the waveforms.
 - If [Auto-save] is enabled, the ECG machine will automatically store the waveforms and reports after printing the waveforms.
- If the waveforms and the reports need to be stored to a specified memory, go to [System Setting] > [Default Memory] and select the desired memory: [Internal Memory] or [USB Flash Drive].

4.3.1 Auto Mode

In Auto Mode, the ECG machine can automatically print the waveforms and reports.



The ECG machine can analyze the resting ECG and output the measurement data, median beat and analysis result, etc.

Simple Report includes patient information, simple measurement data and Minnesota code;

Median Beat Report includes patient information, simple measurement data, Minnesota code, Median Beat waveform and Rhythm waveforms.

In Manual Mode, the analysis reports include analysis report (1), and analysis report (2).

In [Print Setting], when [Analysis Output] is enabled, all the above-mentioned reports will include the analysis results.

After printing waveforms and report, press to print waveforms and report of the last patient.

Auto Print

In Economic mode and Automatic mode, select [Menu] > [Print Setting] > [Auto Print] > [Enable], the ECG machine wil automatically print the waveforms and reports when the following 3 conditions are satisfied.

- Lead off does not detected for 2s (none of any electrodes);
- Five or more QRS complex detected;
- Waveform is stable, no EMG interference or Baseline drift.

EXPLANATION

- The ECG machine analyzes the latest 10 seconds waveforms.
- If the age of the patient is not inputted, the ECG machine will assume this patient as adult during analyzing.

WARNING

- In the case for some special populations (such as pregnant women, the user of vascular drugs, etc.) or mixed by obvious interference in the recording process, the analysis result may be inaccurate. Therefore the final conclusion should be drawn by a doctor, based on analysis result, the clinical characterization of patients and other diagnostic test results.
- If there is too much AC and EMG interference, the identification of P wave and Q wave is not reliable sometimes; if there is baseline wander, the identification of ST segment and T wave is not reliable sometimes.
- If the ending points of S wave and T wave are winding and not clear, it might cause measurement error.
- If R wave is undetected because of low voltage for QRS complex, it might cause some deviations in heart rate measurement.
- If QRS complex has low voltage, the electrical axis measurement and the identification of QRS dividing point can be unreliable.
- Occasionally, the frequent (repetitive) ventricular premature beat might be detected as the median beat.
- When several kinds of arrhythmia occur simultaneously, the identification of P wave might be difficult, and the relative parameters might be unreliable.

4.3.2 Manual Mode

In Manual Mode, press



to start or stop printing. You can switch from one column to another by

pressing



at any time to control the waveform length of every leads.

4.3.3 Upload Mode


In upload Mode, select [Upload] or press to analyze and then upload the ECG waveforms and reports to ECG workstation.

The Printer is disabled in this mode. You can switch to Auto mode or Manual mode to print the reports.

4.4 Advanced Mode

4.4.1 Rhythm Mode

Operate as follows to enter Rhythm Mode:



Enter Rhythm report interface

In the Rhythm interface, the ECG machine begins collecting waveforms of the Rhythm Lead. In single-rhythm pattern, only one lead is selected as rhythm lead and as long as 300s waveform will be collected and analyzed. After waveforms collecting, the ECG machine will automatically analyze the waveforms and enter report interface.

In rhythm report interface, select [Page up] and [Page down] to review more information; select [Print], [Save], [Upload] to print/save/upload the Rhythm report.


Please refer to *ECG Setting* to set [Rhythm Lead] and [Rhythm Time].

EXPLANATION

During the waveform collection, wait for 8 seconds and then press [R-R] to enter report interface.

4.4.2 Analysis Report Mode

Operate as follows to enter analysis report mode:



Enter analysis report

Report (1) includes simple measurements, Minnesota code; median beat waveforms, analysis result and rhythm waveforms.

Report (2) includes waveforms of all the leads.

Analysis report interface is as follows:




Figure 4-3 Analysis Report 1




Figure 4- 4 Analysis Report 2

In above two interfaces, you can select to Upload, Save, and Print the report.

In analysis report (1), **Analysis Result** can be manually edited by the user. Please refer to *List of Interpretation*Codes and Corresponding Description for the details of the analysis result.

Chapter 5 Setting System Parameters

In the Main interface, select $\left[Menu\right]$ to enter the configuration menu.




Figure 5-1 Configration Menu

Press to return to the previous page; press step by step to return to the Main interface.

Select [Save] before exiting to avoid loss of settings caused by sudden power loss.

5.1 ECG Setting

Enter [ECG Setting] to adjust parameters about Electrocardiograph. See following table:

Name	Value	Default	Description	
Lead Mode	Standard Leads	Standard	Cabrera is disabled in iE 101/300.	
Lead Wode	Standard Leads	Leads	Captera is disabled in its 101/300.	
Low page Eilton	25 Hz, 35 Hz, 75 Hz, 100	100 Hz	Calast and setting for large and Eilen	
Low-pass Filter	Hz, 150 Hz, 250 Hz	100 HZ	Select one option for low-pass Filter,	
Baseline Wander	0.01 Hz, 0.02 Hz, 0.05 Hz,	0.5 Hz	Baseline Wander Filter and AC Filter.	

		I	
Filter	0.35 Hz, 0.5 Hz, 0.8 Hz		
AC Filter	OFF, 50 Hz, 60 Hz	50 Hz	
Rhythm Mode	Single-rhythm	Single-rhythm	Three-rhythm is disabled in iE 101/300.
Rhythm Lead	I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6	П	Select a lead as rhythm lead.
Rhythm Time	30 s∼300 s	60 s	Select one option for rhythm recording time.
Wave Sample Time	10 s∼24 s	10s	Select one option for waveform sampling time.
Printout Duration	3~10s	3s	Select one option for printout duration. Note that this option is only available for iE 101.
Pre-acquisition	Disable, Enable	Disable	Set the pre-acquisition mode. If enabled, it can print out pervious waveforms.
Auto-save	Disable, Enable	Enable	Set whether to automatically save the report.
Data Format	ECG	ECG	The data format for storage.
Auto Upload	Disable, Enable	Disable	Set to upload the waveforms and reports automatically or not after printing.
QTC Formula	Bazett, Fridercia, Framingham, Hodges	Bazett	Select one option for QTC formula.
Examination Type	Normal, Physical Examination	Normal Select one option for examination to	

EXPLANATION

For physical examination of a large population, it is recommended to select [Examination Type] > [Physical Examination].

5.2 Print Setting

Enter [Printer Setting] to adjust parameters of the Printer. See following table:

Name	Value	Default	Description	
Gray Level	1~8	6	Select one option for gray level.	
Baseline Width	Baseline Width 1~4		Select one option for waveform thickness.	
Print Speed	5 mm/s, 6.25 mm/s, 10 mm/s, 12.5 mm/s, 25 mm/s, 50 mm/s	25 mm/s	Select one option for paper speed.	
Report Format	Simple Report, Median Beat Report, Manual Report, Disable	Simple Report	Select one option for report format.	
Analysis Output	Enable, Disable	Enable	Select whether to print the interpretation automatically.	
Print Grid Enable, Disable		Disable	Select whether to print the grid on the paper.	
Printer Built-in Thermal Printer, Disable		Built-in Thermal Printer	Select whether to print using built-in printer.	
Record Mode	Refer to the <i>Appendix B</i> Technical Specification and the machine your purchased.	Auto	Select one option for report mode.	
Record Format	Refer to the <i>Appendix B</i> Technical Specification and the machine your purchased.	1 X 12	Select one option for report mode.	
Auto Print	Disable, Enable	Disable	Auto Print is disabled in iE 101/300.	
Print Data Type	Real Time, Synchronous	Real Time	Select one option for print data type.	

5.3 Display Setting

Enter [Display Setting] to adjust parameters about display. See following table:

Name	Value	Default	Description
Display Style	Classic White, Classic Black	Classic White	Set the display style of the screen.
Background	Disable Fredd	E11-	Set to display the background gird or
Grid	Disable, Enable	Enable	not.
		Same Screen	Select to display all the waveforms in
Display Format	Same Screen, Split-screen		one screen or split the waveforms and
			display them in different screens.
Land Farman	Refer to Appendix B	2 V 4	C-14
Lead Format	Technical Specification.	3 X 4	Select one option for lead format.
Lead Standard	Lead Standard IEC Standard, AHA Standard		Select a lead standard to display.

5.4 Patient Information Setting

Enter [Patient Information] to input the patient information.

Input the following information: Sub-ID No., Gender, Age, Date of birth, Height, Weight, BP, Race.

EXPLANATION

Age and Date of Birth cannot be selected at the same time.

5.5 System Setting

Enter [System Setting] to adjust the system parameters. See following table:

Name	Value	Default	Description
			Select normal ECG to display in the
Demo Mode	Normal ECG, Disable	Disable	interface as a demo. Or disable the Demo
			mode.

System Language	English, 中文, etc.	To be determined by the shipping country	Set the system language.
System Version	Version No., Compile Time, Lead Wire Version	/	Display the details of software version.
System Time	Current Time, Date Format, Date, Time	/	Display the details of time and date and to set the date format.
Network Setting	Cable Network	Cable Network	Select one option for network setting.
Transfer Protocol	ТСР	ТСР	The Default setting is TCP. FTP is disabled in iE 101/300.
Cable Network	IP Address, Subnet Mask, Default Gateway	/	Set the value of IP Address, Subnet Mask, and Default Gateway.
Server Setting	IP Address, Port	/	Input the value of IP Address, Port of server.
Silent Mode	ent Mode Disable, Enable		Set to disable or enable the silent mode. If silent mode is enabled, all sound, including alarm tone and key tone will be mute.
QRS Tone			If all three options are zero, a mute icon
Alarm Tone	0~10	6	will be displayed in the Main interface.
Key Tone			
Default Memory Internal Memory, USB Flash Drive		Internal Memory	Select the default memory way for the saved file.
Memory Format Internal Memory Formatting		Internal Memory Formatting	Format the specified memory. The files cannot be recovered after formatting.

Standby Time None, 5 min, 10 min, 30 min, 1 h, 2 h		None	Set the standby time.
Auto Power-Off	None, 30 min, 1 h, 2 h, 3 h	None	Set the auto power off time.
System Password Disable, Enable		Disable	Select to disable or enable to set a system password.
Password Setting 0∼9999		1234	Set the password when the system password is enabled.
General Setting	1~10	1	Select one option, and then set up according to your habit, all your setup will be stored in this option to facilitate your next use. Different doctor or different check up can occupy different options.
Import Setting	Import from USB flash disk	/	Import files from USB flash disk.
Export Setting	Export to USB flash disk	/	Export files to USB flash disk.
Factory Default /		/	Restore to the factory default settings.
Hospital	Hospital /		Input the name of hospital.
Device No. / Input the number of this ECG		Input the number of this ECG machine.	

A CAUTION

Demo mode is designed for representation only. Do not use this mode in clinical analysis, for demo waveforms may be mistaken as that of patient and misdiagnose may happen.

5.6 Factory Maintain

Only the authorized service engineer can enter [Factory Maintain], please contact with our Customer Service Department if necessary.

Chapter 6 Data Management

Select [Menu] > [Data Management] to enter the Data Management interface. Select one source of storage medium from Local (ECG machine), USB 0, or USB 1, the ECG files will be uploaded.




Figure 6- 1 Data Management

No.	Name	Description
1	Select All Shortcut Button	Press to select all the ECG files in the current page.
2	File Attributes Area	Select any file attribute, and sorting files by pressing ENTER .
3	ECG Files Display Area	Display basic patient information of all the ECG files.
4	Open	Open one patient's ECG file.
5	Upload	Upload selected ECG files to ECG workstation or server.
6 and 7	Page Up/Page Down	Browse ECG files in previous or next page.
8	Operate	Select to Refresh, Copy, Move, Delete, or Search the ECG files.

		Search ECG files according to ID number, name, age, time and
		symptom. Search ECG files according to ID number, name, age,
		time and symptom.
		Set one option for the storage medium, including local, USB 0 and
9	Select Storage Medium	USB 1.
10	Indication Area	Indicate pages of ECG files and internal memory.

6.1 Open an ECG File

Select an ECG file, select [Open] and press ENTER to open the file.

EXPLANATION

When you select more than one file to open, the default file is the first selected file.

6.2 Edit an ECG File

After opening the ECG file, you can edit the patient information and analysis result, can also upload, save and print the file. You can refer to *Analysis Report Mode* to know the content and function of the ECG file.

6.3 Delete ECG Files


Select [Operate] > [Delete] and press ENTER to delete files.

≜ CAUTION

Deleted files cannot be recovered. Please confirm before deleting.

6.4 Copy and Move ECG files

Select one or more ECG files, operate as following to transfer the ECG files.



ECG files are copied or moved to the selected medium

EXPLANATION

- Files can be copied or moved between the local ECG machine and external storage medium. ECG files on the local machine will be deleted if the user moves them.
- After selecting files and a path, the files will be copied to the selected path. When copy files to USB device,

the system will create a new folder in the USB device to store the selected files, for example "ecg database".

- When there is not enough memory, the system will indicate that. You should choose new memory to make files copying or moving successfully.
- Please regularly clear data in a storage medium or the speed of the machine will be slowed down.

CAUTION

- When copying or moving the files, the continuity of power supply must be assured, or the files may be missing.
- When copying or moving the files, it is not allowed to insert or unplug the USB device; otherwise, it may cause abnormality of the ECG machine.

Chapter 7 Maintenance

7.1 Main Unit

A CAUTION

- Gently disconnect the Patient Cable and power cable without forcibly pulling the lead wires.
- Clean the ECG machine as well as the accessories and cover the machine from dust.
- Store the machine in a dry and cool environment and avoid excessive shocking and vibration.

7.2 Patient Cable

△ CAUTION

- The lead wire must be periodically checked. Damage may cause incorrect of ECG waveform at some or all leads.
- The user should avoid twisting the Patient Cable, or the life time of the Patient Cable will be shortened.

7.3 Cleaning and Disinfection

Before cleaning, power off the ECG machine and disconnect it from the AC power.

Do cleaning first before disinfection.

The process to clean and disinfect the ECG machine, cables, lead wires and reusable electrodes are as follows:

- Use a clean soft cloth absorbing an amount of cleanser or disinfectant to wipe the surface carefully and avoid touching connectors of the ECG machine and accessories.
- 2) When necessary, wipe the superfluous cleanser or disinfectant with dry cloth.
- 3) Place in the ventilated and cool environment to dry the ECG machine and accessories.

Sterilization operation for this ECG machine and accessories is not recommended, unless the manual of the accessories has requirement.

△ CAUTION

■ While cleaning and disinfection, do not splash liquid into the ECG machine and the accessories.

- Disinfections may cause damage to the ECG machine or accessories to a certain degree. It is suggested that only when necessary, disinfect the ECG machine and accessories.
- Neutral cleanser or disinfectant is recommended.

7.4 Recording Paper

△ CAUTION

- To ensure good ECG recording, please use suitable thermal recording paper for the ECG machine. Incorrect recording paper can damage printer head and cause problems such as blurring trace and incorrect paper running. Pay attention to the following comments on recording paper.
- Never use recording paper coated with wax for the ECG machine. It may cause serious problems to the printer head.
- When exposed to high temperature, high humidity and direct sunlight, the recording paper will deteriorate. It is therefore required to store the thermal recording paper in a dry and cool environment.
- When exposed to fluorescent light for long time, the recording paper will deteriorate.
- When stored with polyvinyl chloride (PVC), the recording paper will deteriorate.
- If the thermal recording paper is stored overlapping for a long time, the printing impression will leave traces in other pages, which will cause mislead readings.
- Use suitable size recording paper for the ECG machine. Or, it may cause damage to the printer head and Silicon rubber shaft.

7.5 Battery

The ECG machine is equipped with a built-in rechargeable battery to assure continuous operation when AC power is unavailable. Charging, capacity indication and replacement of the battery are described below:

♦ Charging

The ECG machine is designed with a charger and protector for the battery.

- Please turn off the machine first before charging the battery.
- The battery-charging indicator on the operation keyboard will become green when the battery is charged completely.

■ Discharge and charge the battery at least once every three months (discharge the battery until the machine turns off automatically, and fully charge the battery).

♦ Capacity indicator

When the unit is powered by battery, there will be a symbol of battery capacity indication displayed on the LCD.

For example:



Full battery capacity, it can work continuously for about 3 hours.



Battery capacity is sufficient.



Insufficient battery capacity, charging is required.



Battery capacity is going to running out, immediate charging is demanded.



Battery capacity has already run out and blackout may happen at once, immediate charging is demanded.

♦ Battery replacement

The battery should be replaced by the professionals according to the following procedures.

- 1. Power off the ECG machine and disconnect the AC power cable.
- 2. Flip over the ECG machine and disassemble the battery back cover based on the instruction on the back cover.
- 3. Disconnect the battery plug and take out the battery.
- 4. Replace the existing battery with a new one. Pay attention to polarity and connection.
- 5. Install the battery cover.




Figure 7-1 Battery Replacement

EXPLANATION

Refer to *Foreword* for other warning information of battery.

7.6 Silicon Rubber Shaft for Printing

The silicon rubber shaft shall be kept clean, smooth and free from dirt. Otherwise, the ECG machine may print out unsatisfied ECG trace. To clean the shaft, wipe the shaft with soft cotton moistened with medical alcohol and at the same time rotate shaft until it is clear enough.

7.7 Thermal Print Head

Residue and dirt on the thermal print head could affect the clarity of recorded ECG waveform. To clean the thermal print head, open the paper cover and clean the print head with soft cotton moistened with medial alcohol. It is not permitted to clean the print head with a sharp object, which can cause permanent damage to the print head. Thermal print head maintenance should be done at least once a month.

Chapter 8 Troubleshooting

8.1 Lead Fault

1. Data saturation or overflow happens.

Solution:

Ensure that all leads are in good contact, and wait for half a minute or the waveforms on the screen are stable, and then start printing.

2. Straight line is printed in some leads.

Solution:

- Check if the metal piece of limb electrode contacts the body properly; if not, adjust the position of the limb electrode, and adjust the tightness if necessary.
- 2) Check if the limb electrodes and chest electrodes are oxidized or faded, and clean the accessories or replace with new ones. Oxidation and aging cause conductive deterioration of the electrode sheet, resulting in poor signal transmission.
- 3) Treat the skin of the patient with alcohol; especially in the winter of the North, dry skin causes skin resistance to become larger, which will impact the signal reception.
- 4) Please clean the joints of lead wires, suction ball and limb clip, reinstall and tighten all joints. After long-term use, the joints will have dirt or become loose, resulting in poor signal transmission.
- 5) Check if the appearance of the lead wires has obvious fracture; if yes, replace with new lead wire. If not, connect a proper lead wire to the device. If the waveform is stable, the lead wires have problem and have to be replaced.
- If there is no lead wire available, check if the lead wires conduct with a multimeter. First check if the inner conductors of the lead wires are conductive. Generally speaking, the acceptable resistance shall be about $10 \text{ k}\Omega$. Then check if there is a short circuit between the outer shield and inner conductor. The resistance shall be infinity. If the lead wires have a problem, please contact our service department to replace new lead wires.
- The other reason for lead fault can be caused by failure in signal communication. Please exclude other causes for lead fault problems first, and contact our service technicians if necessary.

8.2 Printer Failure

1. Unclear printing.

Solution:

- 1) Whenever a printer fault occurs, such as poor or incorrect ECG recording, you may try to clean the thermal printer head with soft cotton dipped with medical alcohol.
- 2) Please check if the quality of the thermal paper is poor or if the unsealing time is too long, resulting in reduced performance of thermal layers, and replace with provided or specified recording paper.
- 3) If the above methods are not applicable, guide the user to test the print head and check if the print head has breaking point; if yes, contact the company service department to replace the thermal print head.
- 2. Upper half or lower half is blank.

Solution:

Check if the bearing on both ends of the rubber shaft of the paper compartment cover is worn, and replace with new bearing if yes.

3. All paper or most part is blank.

Solution:

- 1) Make sure that the thermal recording paper is not installed backwards.
- Check if the print head is stuck by dirt (such as adhesive tape); this often occurs when new print paper is replaced.

8.3 Indicator of Lead Off

This ECG machine can detect the lead connecting status continuously. When the leads are not well connected to the main unit, it means that the signals cannot be transferred correctly, thus there is "lead off:*" indication, accompanied by voice alarm. The symbol "*" represents the fault lead, the waveform of which will display as a straight line. Please check carefully whether the connection among the related electrodes, human body, lead wire and the main unit remains well.

8.4 AC Interference

Apparent and regular trembling of ECG waveform in the process of recording due to AC interference is shown as below.



Causes of baseline wander are varied, please do following checks one by one:


- 1. Make sure the ECG machine is properly grounded according to the instructions.
- 2. Make sure the lead wires or electrodes are properly connected.
- 3. Make sure the electrodes and the patient skin have been covered with ECG gel.
- 4. Make sure the exam bed is properly grounded.
- 5. Make sure the patient is not in touch with the wall or the metal part of the bed.
- 6. Make sure the patient is not in touch with anybody else.
- 7. There shall be no large power electric equipment (such as x-ray machine, ultrasound scanner and so on) operating nearby.
- 8. The patient shall not be wearing some jewelry like diamond.

A CAUTION

Set AC filter to ON if AC interference still exists after the above checks are completed.

8.5 EMG Interference

Irregular trembling of ECG waveform due to EMG interference is shown as below.



Causes of baseline wander are varied; please do following checks one by one:

- 1. Make sure that the exam room is comfortable for examination.
- 2. Soothe the patient from irritation or excitement.
- 3. Make sure the exam bed is in suitable size.
- 4. Never have talks with the patient during ECG recording.


5. Make sure the limb electrode is too tight to make the patient uncomfortable.

A CAUTION

Set EMG filter to lower value if EMG interference still exists after the above checks are completed, and keep in mind that the recorded ECG waveform, particularly R wave, will be somewhat attenuated.

8.6 Baseline Wander

Irregular movement of ECG baseline due to baseline drift is shown as below:



Causes of baseline wander are varied, please do following checks one by one:

- 1. Make sure the electrodes are in good contact with skin.
- 2. Make sure the connection between the Patient Cable and electrodes is normal.
- 3. Make sure the electrodes are clean and patient skin contacted with the electrodes has been well pretreated.
- 4. Make sure the electrodes and skin are covered with ECG gel.
- 5. Keep the patient silent and motionless, and keep the patient from hyperventilation.
- 6. Used electrodes shall not be utilized with new ones in patient examination.

A CAUTION

If the problem still can't be cleared, please turn up the baseline wander filter, keep in mind that the recorded ECG waveforms, particularly T wave and ST segment, will be somewhat distorted.

8.7 The ECG Machine cannot be turned on

1. AC power is not working properly and the battery is exhausted.

Solution:

First check if the power outlet connects properly, if the power line and the machine connect properly, and then check if local AC voltage is normal. If everything is ok, check if the fuse is good. If everything is

-- 54 -- User Manual for Electrocardiograph

normal after above examination, return the machine to the manufacturer for repair, for it may be damaged.

2. After turning on, the machine turns off automatically after a few minutes.

Solution:

While working, if the screen displays the battery power, the machine is using the battery, but the battery power is insufficient, resulting in automatic shutdown. Please supply the machine with AC power, or charge the battery before working; if the user is using AC power and the machine still turns off automatically, please check according to step1.

8.8 Paper Feeding Failure

1. Press but paper is not fed.

Solution:

Check if the keyboard has been damaged and replace damaged keyboard.

2. Press , the paper is not fed, and there is abnormal sound from recorder.

Solution:

First check if the recording paper is installed properly, and if the gear on the cover of paper magazine is in good condition.

3. Paper feeding isn't smooth, paper is stuck, or waveforms are compressed.

Solution:

First check if the thermal paper complies with the standard, then if the paper is installed properly, or replace a new roll of paper. Finally, replace the paper shaft.

8.9 Battery is quickly Charged and Discharged

If the battery is often not fully charged, the performance will be deteriorated.

Solution:

It is recommended to charge the battery continuously until the battery is fully charged and activated for the first two times. Supply the machine with AC power as far as possible.

8.10 Wrong Analysis Result

For the case for some special populations (such as pregnant women, the user of vascular drugs, etc.) or mixed by obvious interference in the recording process, the analysis result of the resting ECG analysis of this ECG machine

User Manual for Electrocardiograph -- 55 --

may be inaccurate. The possible reason may be as followings:

- 1. Poor contact between electrode and patient skin, caused by improper skin treatment and incorrect connection.
- 2. The patient has relatively large movement in the recording process.
- 3. Gender or age isn't entered;
- 4. If there is too much AC, EMG and breathing interference, the identification of P wave and Q wave is not reliable sometimes; if there is baseline wander, the identification of ST segment and T wave is not reliable sometimes.
- 5. If QRS complex has low voltage, R-wave may be missed, and the electrical axis measurement and the identification of QRS dividing point can be unreliable. Or frequent ventricular contraction occurs or a variety of arrhythmias merge, the relevant detection parameters may be unreliable.
- 6. The filter settings are incorrect.

Solution:

- 1) Treat as *Apply Electrode* and wait until the waveform is stable before reanalyzing.
- 2) Enter the patient gender and age correctly.
- 3) Exclude the interference as the methods described in AC Interference, EMG Interference and Baseline Wander before reanalyzing.
- 4) Reset to an appropriate filtering value.

8.11 File Uploading Failure

The most possible reason is that network settings have problems, please check the network connection and refer to *Network Connection* to re-set the network.

Appendix A Package and Accessories

A.1 Packing List

Туре	Item	Qty
	ECG Machine	1 unit
	Lead Wires	1 set
	Limb Electrodes	1 set
	Chest Electrodes	1 set
	Power Cable	1 piece
G. 1 1	Equipotential cable	1 piece
Standard	Thermal Recording Paper	1 piece
	Paper Roller	1 piece
	User Manual	1 copy
	Warranty Card	1 сору
	Qualified Certificate	1 copy
	Packing List	1 сору
	Electrodes	10 piece
	Pediatric Limb Electrodes	1 set
Optional	Pediatric Chest Electrodes	1 set
	Clips	1 piece
	Cardiograph Cart	1 unit

A.2 Dimensions and Weight

	281 mm X 191 mm X 59 mm (iE 101)
Length × Whith × Height	281 mm X 191 mm X 59 mm (iE 300)
Not Weight	About 1.3 kg (iE 101)
Net Weight	About 1.3 kg (iE 300)

A CAUTION

- Please open the package according to instructions on the packing box.
- Enclosed accessories and documents shall be checked according to the packing list before starting checking on the unit.
- Whenever there will be mismatch of the accompanying materials with the packing list, contact our Customer Service Department immediately.
- To ensure good performance and safe operation of the ECG machine, please use the accessories supplied by the manufacturer.
- The package box should be kept well for the regular inspection or maintenance for the machine.

Appendix B Technical Specification

B.1 Specifications

B.1.1 Main Unit

Lead	Standard 12-lead
Acquisition Mode	Simultaneous 12-lead Acquisition
D 15	1 X 12, 1 X 12+1R (iE 101)
Record Format	1 X 12, 1 X 12+1R, 3 X 4, 3 X 4+1R, 3/2 (iE 300)
Record Mode	Auto, Manual, Upload
Display Format	3 X 4, 3 X 4+1R, 6 X 2, 6 X 2+1R, 12 X 1
Rhythm Time	30~300s waveforms acquisition for rhythm analysis
Measurement	Ventricular Rate, PR Interval, QRS Time Limit, QT/QTC Interval, P/QRS/T Axis,
Parameters	RV5/SV1 Amplitude and RV5+SV1 Amplitude
	AC Filter
T	Baseline Wander Filter
Filters	EMG Filter
	High-frequency filter
Input CIR Current	≤0.1 µA
Input Impedance	≥30 MΩ (Full-band)
Time Constant	≥3.2 s
Frequency Response	0.05 Hz∼250 Hz
Noise Level	\leq 12.5 μ V _{p-p}
Sensitivity Threshold	\leq 20 μ V _{p-p}
Sensitivity	Auto, 0.625 mm/mV, 1.25 mm/mV, 2.5 mm/mV, 5 mm/mV, 10/5 mm/mV, 10
Sensitivity	mm/mV, 20/10mm/mV, 20 mm/mV, and 40 mm/mV
Standard Sensitivity	$10 \text{ mm/mV} \pm 2\%$
Calibration Voltage	1 mV±3 %

	Using the method described in 4.2.7.1 of AAMI EC11 to test the overall system error,	
	which is within ±5%;	
	Using method A and D described in 4.2.7.1 of AAMI EC11 to test frequency response.	
Accuracy of input	Because of sampling characteristics and the asynchronism between sample rate and	
signal reproduction	signal rate of the ECG machine, digital systems may produce a noticeable modulating	
	effect from one cycle to the next, particularly in pediatric recordings. This	
	phenomenon, which is not physiologic, shall be clearly described in the operator's and	
	service manuals.	
CMRR	>115 dB	
Patient Leak Current	<10 μΑ	

B.1.2 Recorder Specification

Recorder	Thermal Dot Matrix Word Printing System	
Danada a Danas	50 mm, roll paper (iE 101)	
Recording Paper	80 mm, roll paper (iE 300)	
Paper Speed	(5, 6.25, 10, 12.5, 25, 50)mm/s ± 5%	

B.1.3 Other Specification

Display on LCD	5-inch TFT LCD screen (iE 101)	
Display on LCD	5-inch TFT LCD screen (iE 300)	
Safety Classification	IEC60601-1, Class I, Type CF	
A C D C I	100 V-240V, 50 Hz /60 Hz, 80 VA (iE 101)	
AC Power Supply	100 V-240V, 50 Hz /60 Hz, 80 VA (iE 300)	
	Rechargeable lithium battery, 11.1 V/ 2600mAh.	
	In environment temperature 25 $^{\circ}$ C \pm 5 $^{\circ}$ C and with the machine turning off, the	
DC Power Supply	charging time is not more than 2 hours to charge the battery to 90%.	
	In environment temperature 25 $^{\circ}$ C \pm 5 $^{\circ}$ C, the continuous working time is not less than	
	3 hours while the ECG device is continuously printing.	

B.2 Environment Requirements

	Transportation				
1	Environment Temperature	-20 °C∼+55 °C			
	Relative Humidity	≤95 % (No condensation)			
	Air Pressure	70 kPa∼106 kPa			
	Transportation: avoid direct sunshine and rain.				
	Storage				
2	Environment Temperature	-20 °C∼+55 °C			
	Relative Humidity	≤95 % (No condensation)			
	Air Pressure	70 kPa∼106 kPa			
	The packed ECG should be stored in the well-ventilated room without corrosive gases.				
	Using				
3	Environment temperature	+5 °C ∼+40 °C			
	Relative humidity	≤95 % (No condensation)			
	Air pressure	70 kPa∼106 kPa			

-- Blank Page --

Appendix C Working Principle and Block Diagram

C.1 Power Supply Subsystem


1) Working principle:

By ways of AC to DC high frequency power switching technique, the switching power supply output 20V DC voltage. This voltage works to charge the built-in rechargeable battery at constant voltage and limited current, and at the same time, applies to the power switchover circuitry together with the battery output. If the switchover circuitry is turned on, several stable power outputs will be generated through various switching power technique, including main power supply of +5V and +24 V (+12 V) by the switching power stabilizer, +3.3 V, +1.8 V and -5 V by the power supply transformation.

+5V, +3.3 V, +1.8 V and -5 V supply power for the control system circuitry with the equipotential as reference with 750mA normal load, 3A output current capacity and 3.75A output current limit for short circuit protection.

The +24 V (+12 V) voltage supplies power for the paper driving motor and the thermal printer. The motor is driven by means of width modulation and wave chopping technique in sequence to improve power efficiency. It has about 500mA as normal load, 850mA as output current capacity and 1.2A output current limit for short circuit protection. The self-excitation power switching circuitry transforms the output of the switchover circuitry to several voltages for analog circuitry. The +5V voltage needed by isolated digital circuitry is the direct stabilized output of the switching power supply, of which the normal load is about 150mA and the current capacity is 300mA. The +8V and -8V voltage for the isolated analog circuitry is the un-stabilized output from the switching power supply. Its normal load is about 60mA and has 100mA current capacity.

2) Block diagram (Schematics and parts of list of this unit are only provided for qualified service center under supervision of the manufacturer.)




C.2 Acquisition Module

1) Working principle

Acquisition Module is connected to the main unit by ECG port which provides power supply and communication interface. When connecting, the analog circuit and control circuit of acquisition module are powered through the isolated DC/DC switching output. The analog circuit is composed of input protection circuit, anti-aliasing low-pass filter and ECG analog front chips. The electrodes acquire the millivolt electrical signals from human body, which will be converted to digital signals by ECG analog front chips first, then transferred to the processor that controls the converter in the ECG analog circuit, ECG data processing and operation keys on acquisition module. The sampling rate of signals is 8000Hz with 250Hz as bandwidth (-3dB), which meets the standard of AHA and CSE (sampling rate no less than 500Hz). After the signals are processed and filtered, they are sent back to the main unit through optically coupled isolation interface.

2) Block diagram:



3) Leads of acquisition module:

Lead nomenclature	Definition	Name of lead
I	I =L-R	Diamalan limb laada
II	II =F-R	Biopolar limb leads
III	III =F-L	(Einthoven)
aVR	aVR=R- (L+F)/2	Assembled looks
aVL	aVL=L- (R+F)/2	Augmented leads
aVF	aVF=F- (L+R)/2	(Goldberger)
V1	V1=C1- (L+R+F)/3	
V2	V2=C2- (L+R+F)/3	
V3	V3=C3- (L+R+F)/3	Unpolar chest leads
V4	V4=C4- (L+R+F)/3	(Wilson)
V5	V5=C5- (L+R+F)/3	
V6	V6=C6- (L+R+F)/3	

C.3 Control System


1) Working principle

The control system consists of controllers for printer, keyboard, LCD screen and a CPU subsystem. Through high-speed photo-couplers, CPU subsystem receives lead signals from the Data acquisition Subsystem and applies them to the printer controller after digital filtration, gain adjustment and printers driving then

User Manual for Electrocardiograph -- 65 --

complete the ECG waveforms printing. The lead data will also be measured and interpreted by CPU. In addition to measurement and interpretation on the printed ECG waveform, CPU also receives interruption signals and key codes from the keyboard controller to carry out key interrupt routine. Acquisition and processing of signals for detection of leadoff, out-of-paper detection, battery capacity management, automatic power off, CRO output and EXT input are all managed by CPU. The printer controller receives instructions and data from the CPU, and works to manage the buffering area and generate control signals for the stepping motor and thermal print head to print out ECG waveforms and related information. The keyboard controller works to generate keyboard scan signals, removes key bounces when key is pressed, and sends key codes and interruption signals to CPU for further processing. LCD controller receives instruction s and data from CPU, and works to display the unit's operation status

2) Block diagram



Appendix D List of Interpretation Codes and Corresponding Description

Arrhythmia
Marked rhythm irregularity
Sinus rhythm
Sinus arrhythmia
Marked sinus arrhythmia
Sinus tachycardia
Sinus bradycardia
Atrial rhythm
Atrial fibrillation
Atrial fibrillation with rapid ventricular response
Atrial fibrillation with slow ventricular response
Atrial fibrillation with aberrant conduction, or ventricular premature complexes
Atrial fibrillation with rapid ventricular response with aberrant conduction,or ventricular
premature complexes
Atrial fibrillation with slow ventricular response with ventricular premature complexes
Atrial tachycardia
Atrial flutter
Atrial flutter with aberrant conduction or ventricular premature complexes
Cannot rule out atrial flutter
Junctional rhythm
Junctional tachycardia
Supraventricular rhythm
Supraventricular tachycardia
Supraventricular bradycardia
with occasional supraventricular premature complexes
with frequent supraventricular premature complexes
with frequent supraventricular premature complexes in a pattern of bigeminy
Ventricular rhythm

Ventricular tachycardia with occasional ventricular premature complexes with frequent ventricular premature complexes with frequent ventricular premature complexes in a pattern of Bigeminy Undetermined regular rhythm Undetermined rhythm with occasional ectopic premature complexes with frequent ectopic premature complexes with frequent ectopic premature complexes in a pattern of Bigeminy Electronic Atrial pacemaker Electronic ventricular pacemaker Dual chamber Electronic pacemaker Demand pacemaker Pacemaker failure Ventricular Fibrillation **AV Conductive Defect** Possible third degree AV block Second degree AV block, Wenckebach type Second degree AV block, Mobitz type II First degree AV block Short PR interval Type-A WPW syndrome Type-B WPW syndrome Atypical WPW syndrome Intermittent WPW syndrome **Intraventricular Conductive** Left bundle branch block Incomplete left bundle branch block Right bundle branch block, plus RVH Right bundle branch block Incomplete right bundle branch block

RSR' in lead V1/V2, consistent with right ventricular conduction delay Left anterior fascicular block Left posterior fascicular block Nonspecific intraventricular conduction block Nonspecific intraventricular conduction delay Myocardial Cannot rule out anterior myocardial infarction, probably old Cannot rule out anterior myocardial infarction, age undetermined Possible anterior myocardial infarction, possible acute Possible anterior myocardial infarction, probably recent Possible anterior myocardial infarction, probably old Possible anterior myocardial infarction, age undetermined Anterior myocardial infarction, possible acute Anterior myocardial infarction, probably recent Anterior myocardial infarction, probably old Anterior myocardial infarction, age undetermined Cannot rule out anteroseptal myocardial infarction, probably old Cannot rule out anteroseptal myocardial infarction, age undetermined Possible anteroseptal myocardial infarction, possible acute Possible anteroseptal myocardial infarction, probably recent Possible anteroseptal myocardial infarction, probably old Possible anteroseptal myocardial infarction, age undetermined Anteroseptal myocardial infarction, possible acute Anteroseptal myocardial infarction, probably recent Anteroseptal myocardial infarction, probably old Anteroseptal myocardial infarction, age undetermined Cannot rule out anterolateral myocardial infarction, probably old Cannot rule out anterolateral myocardial infarction, age undetermined Possible anterolateral myocardial infarction, possible acute Possible anterolateral myocardial infarction, probably recent Possible anterolateral myocardial infarction, probably old

Possible anterolateral myocardial infarction, age undetermined
Anterolateral myocardial infarction, possible acute
Anterolateral myocardial infarction, probably recent
Anterolateral myocardial infarction, probably old
Anterolateral myocardial infarction, age undetermined
Cannot rule out septal myocardial infarction, probably old
Cannot rule out septal myocardial infarction, age undetermined
Possible septal myocardial infarction, possible acute
Possible septal myocardial infarction, probably recent
Possible septal myocardial infarction, probably old
Possible septal myocardial infarction, age undetermined
Septal myocardial infarction, possible acute
Septal myocardial infarction, probably recent
Septal myocardial infarction, probably old
Septal myocardial infarction, age undetermined
Cannot rule out lateral myocardial infarction, probably old
Cannot rule out lateral myocardial infarction, age undetermined
Possible lateral myocardial infarction, possible acute
Possible lateral myocardial infarction, probably recent
Possible lateral myocardial infarction, probably old
Possible lateral myocardial infarction, age undetermined
Lateral myocardial infarction, possible acute
Lateral myocardial infarction, probably recent
Lateral myocardial infarction, probably old
Lateral myocardial infarction, age undetermined
Cannot rule out inferior myocardial infarction, probably old
Cannot rule out inferior myocardial infarction, age undetermined
Possible inferior myocardial infarction, possible acute
Possible inferior myocardial infarction, probably recent
Possible inferior myocardial infarction, probably old
Possible inferior myocardial infarction, age undetermined

Inferior myocardial infarction, possible acute Inferior myocardial infarction, probably recent Inferior myocardial infarction, probably old Inferior myocardial infarction, age undetermined Cannot rule out inferior myocardial infarction with posterior extension, probably old Cannot rule out inferior myocardial infarction with posterior extension, age undetermined Possible inferior myocardial infarction with posterior extension, possible acute Possible inferior myocardial infarction with posterior extension, probably recent Possible inferior myocardial infarction with posterior extension, probably old Possible inferior myocardial infarction with posterior extension, age Undetermined Inferior myocardial infarction with posterior extension, possible acute Inferior myocardial infarction with posterior extension, probably recent Inferior myocardial infarction with posterior extension, probably old Inferior myocardial infarction with posterior extension, age undetermined posterior myocardial infarction, possible acute posterior myocardial infarction, age undetermined Abnormal Q wave ? [Lat., Inf.] Abnormal Q wave ? [Ant.] Abnormal Q wave ? [Ant., Lat.] Abnormal Q wave ? [Ant., Inf.] Abnormal Q wave ? **Ventricular Hypertrophy and Atrial** Possible right ventricular hypertrophy Right ventricular hypertrophy Right ventricular hypertrophy, probably repolarization abnormality Minimal voltage criteria for LVH Possible left ventricular hypertrophy left ventricular hypertrophy Left ventricular high voltage (moderate) Left ventricular hypertrophy, probably repolarization abnormality

Possible left atrial enlargement Left atrial enlargement Possible right atrial enlargement Right atrial enlargement Biventricular hypertrophy Biventricular hypertrophy with repolarization abnormality **Axis Deviation** Moderate left axis deviation Abnormal left axis deviation S1-S2-S3 pattern Abnormal right axis deviation Moderate Right axis deviation Indeterminate axis **ST-T Abnormality** ST depression, possible digitalis effect Minimal ST depression Moderate ST depression Marked ST depression, possible subendocardial injury Marked ST depression, possible subendocardial injury or Digitals effect Marked ST depression, consistent with subendocardial injury Junctional ST depression, probably normal Abnormal junctional ST depression Possible anterior injury or acute infarct Anterior injury or acute infarct Possible anteroseptal injury or acute infarct Anteroseptal injury or acute infarct Possible anterolateral subepicardial injury Anteroseptal subepicardial injury Possible septal subepicardial injury Septal subepicardial injury Possible lateral subepicardial injury

Lateral subepicardial injury
Possible inferio injury or acute infarct
Inferio injury or acute infarct
T wave abnormality, possible anterior ischemia
T wave abnormality, possible anterior ischemia or digitalis effect
T wave abnormality, consistent with anterior ischemia
T wave abnormality, possible anterolateral ischemia
T wave abnormality, possible anterolateral ischemia or digitalis effect
T wave abnormality, consistent with anterolateral ischemia
T wave abnormality, possible lateral ischemia
T wave abnormality, possible lateral ischemia or digitalis effect
T wave abnormality, consistent with lateral ischemia
T wave abnormality, possible inferio ischemia
T wave abnormality, possible inferio ischemia or digitalis effect
T wave abnormality, consistent with inferio ischemia
ST elevation, probably early repolarization
Early repolarization
ST elevation, consistent with subepicardial injury, pericardiatis, or Early repolarization
Possible acute percarditis
acute percarditis
Nonspecific ST&T wave abnormality
Nonspecific ST&T wave abnormality, probably digistalis effect
Tall T waves, possible hyperkalemia
Nonspecific T wave abnormality
Nonspecific T wave abnormality, probably digitalis effect
Nonspecific ST elevation
Others
Low voltage
Low voltage in limb leads
Low voltage in chest leads
Long QT interval

Short QT interval

Dextrocardia?

LIMB LEADS REVERSED

Abnormal QRS-T angle

Consistent with pulmonary disease

Artifacts present

Cannot be analyzed, re-record recommended

Overall Judgment

Normal ECG

Borderline ECG

Atipical ECG

Abnormal rhythm ECG

Abnormal ECG


A CAUTION

List of Interpretation codes may be subject to changes without notice.

Appendix E Measurement, Diagnosis, Analysis and Assessment of ECG Machine

E.1 Methods to determine the amplitude of P, QRS, ST and T wave

1) P wave amplitude




P On is the starting position of P wave, P Off is the ending position of P wave, and dashed line is the reference baseline

To measure P wave amplitude: the line from the starting point to the ending point of P wave is the reference baseline, as shown in above Figure. The positive amplitude is from the reference baseline to top edge of P wave; the negative amplitude is from the reference baseline to bottom edge of P wave.

2) QRS complex, ST segment and T wave amplitude


When measuring QRS complex, ST segment and T wave amplitude, the horizontal line of QRS complex beginning part is used as the reference baseline, as shown below:



The measurement of QRS complex, ST segment and T wave amplitude uses the horizontal line of QRS complex

beginning part as the reference baseline QRS Onset is the starting position of QRS wave

E.2 Processing method of isoelectric segment in QRS complex



Isoelectric segment between dash lines are in QRS complex

As shown above, the isoelectric segment beginning from the starting position of QRS complex is processed as a part of QRS complex, but doesn't belong to the meaningful wave later (waveform area is larger than $160 \,\mu\text{V} \cdot ms$)

E.3 Low incidence heart disease not included in testing and diagnosis database

Test with CSE database, but this database doesn't have sufficient number of acute myocardial infarction and myocardial ischemia ECG.

E.4 ECG diagnosis categories and the number of ECG test of each category

The accuracy of disease diagnosis and non-ECG means used to verify the effectiveness of heart disease diagnosis, as well as the patients statistics data (e.g. age, gender, race) of each group.

Test with CSE database, the following table lists disease diagnostic categories, the number of ECG testing of each category and the accuracy of disease diagnosis.

CSE database sample properties are as follows:

Total number of samples: 1220 (male: 831, female: 389)

Race: White

Age: 52 ± 13

Type of disease	ECG test	Sensitivity	Specificity	Positive
	number	(%)	(%)	predictive value
				(%)
Normal	382	92.7	73.9	61.8
Left ventricular hypertrophy	183	60.1	97.0	77.7
Right ventricular hypertrophy	55	32.7	99.9	92.3
Biventricular hypertrophy	53	26.4	99.9	93.3
Anterior wall myocardial infarction	170	80.6	97.7	85.1
Inferior wall myocardial infarction	273	67.0	97.8	89.7
Composite myocardial infarction	73	64.7	99.7	94.0
Hypertrophy and myocardial infarction	31	46.8	100.0	100.0

E.5 The smallest waveform identified by the device and the stability of measurement when noise exists

If the area of certain waveform is greater than or equal to $160\,\mu\text{V}\cdot\text{ms}$, it is considered as meaningful wave, otherwise it is meaningless. Recognizing meaningful waveforms in area method can effectively reduce the noise. The stability of the measurement when noise exists is shown below

Overall measurement parameter	Type of added noise	Mean difference (ms)	Variance (ms)
	High frequency	-0.1	0.64
P time limit	Power frequency	0.25	1.5
	Low frequency	-2.3	3.8
	High frequency	1.6	2.4
PR interval	Power frequency	-0.1	1.5
	Low frequency	0.38	9.5
	High frequency	0.75	4.0
QRS time limit	Power frequency	-1.1	1.7
	Low frequency	0.3	4.4
	High frequency	-1.6	3.6
QT interval	Power frequency	-0.5	1.2
	Low frequency	4.9	5.6

E.6 Low incidence cardiac rhythm not included in the ECG rhythm test database

The low incidence cardiac rhythms not included in the test database:

- 1. Grade II conduction block;
- 2. Grade III conduction block.

E.7 ECG rhythm diagnosis categories and ECG test number of each category

Accuracy of rhythm diagnosis and the patient statistics data (e.g. age, gender, race) of each group

The following able gives the rhythm categories, ECG test number of each category and accuracy of disease diagnosis.

The test database sample properties are as follows:

Total number of sample: 4500 (male: 2847, female: 1653)

Race: Yellow

Age: 48 ± 12

-- 78 -- User Manual for Electrocardiograph

Rhythm type	ECG test	Sensitivity	Specificity	Positive predictive
	number	(%)	(%)	value (%)
Sinus rhythm	3656	98.0	91.1	97.9
Premature ventricular contractions	351	87.2	98.9	81.2
Supraventricular premature beats	247	68.8	99.6	89.9
Atrial fibrillation	192	89.6	98.7	91.0
Atrial flutter	49	65.3	99.9	88.9
Pacemaker rhythm	5	100.0	100.0	100.0

E.8 Sensitivity regularly test instructions

Inspect ECGs: EGC-1C

Inspection methods:

- Make ECG machine set in lead I, the sensitivity is set as 10 mm/mV., EGC-1C transmits the U_{in} as 1 mV, frequency 10 Hz sine wave signal to the ECG machine.
- Test the waveform amplitude h_m on the Inspected ECG machine. Calculate the corresponding to deviations of the sensitivity according to the following formula, should meet the maximum allowable relative deviation of \pm 5 %.

$$\delta_S = \frac{S_{\rm m} - S_{\rm n}}{S_{\rm n}} \times 100\%$$

The formula: S_n - nominal value of Sensitivity;

S_m-test value of sensitivity;

h_m-the waveform amplitude of sensitivity ;

U_{in}-input signal amplitude if the inspected ECG machine

- Make ECG machine set in lead I, the sensitivity is set as 20 mm/mV. EGC-1C transmits the Uin as 0.5 mV, frequency 10 Hz sine wave signal to the ECG machine. Using the same method to test the relative deviation of 20 mm/mV sensitivity.
- 4) Make ECG machine set in lead I, the sensitivity is set as 5 mm/mV. EGC-1C transmits the Uin as 2 mV, frequency 10 Hz sine wave signal to the ECG machine. Using the same method to test the relative deviation of 5 mm/mV sensitivity.
- 5) Make ECG machine set in lead I, the sensitivity is set as 2.5mm/mV. EGC-1C transmits the Uin as 4 mV, frequency 10 Hz sine wave signal to the ECG machine. Using the same method to test the relative deviation

 User Manual for Electrocardiograph -- 79 --

- of 2.5 mm/mV sensitivity.
- 6) According to the 1 and 2 steps to change the leads of the ECG machine, and make the ECG-1C's output signals connected to corresponding lead of the ECG machine, to complete all channel's inspect, and then select the largest relative deviation from the test results for each test point, as the result of the inspection.

E.9 Distortion test

The function of ECG machine will not be affected adversely by the running of the pacemaker, which can be verified in the following way:

- a) Superimpose the pulse wave of 200 mV peak, rise time less than 100 µs, 1ms pulse width, and 100 beats / min repetition rate with the sine wave signal of 1mV peak-valley value and 40 Hz frequency, and input to the ECG machine (set to standard sensitivity). The time required to restore the sine wave signals recorded by the ECG machine to 70 % of the initial value (when peak-valley value is 1mV and gain is 10 mm/mV, the initial value should be 10 mm) shouldn't exceed 50 ms; in the above test, the maximum baseline drift accumulated in 10 s doesn't exceed 10 mm; both with and without pulse, the amplitude difference recorded by sine wave signals (after waveform is stable) isn't greater than ±1 mm.
- b) The filter of ECG machine must be opened for distortion tests.
- c) The ECG machine can pass one of the following two tests:
 - String the pacemaker pulse wave of 200 mV peak, rise time less than 100 μs, 1ms pulse width, and 120 pulses / min repetition rate together with the symmetrical triangular wave of 2 mV amplitude and 100 ms duration. The starting time of pulse wave should be 40 ms earlier (or later) than the starting time of triangular wave, input such a signal to the ECG machine, record in the standard sensitivity, the triangular wave is clearly visible on the ECG machine records, the difference between recorded amplitude and the original amplitude (the original amplitude of the waveform with 2 mV amplitude should be 20 mm under 10 mm/mV gain) does not exceed 20 %, and the location of the pacemaker pulse can be clearly identified in the ECG machine records.
 - String the pacemaker pulse wave of 200 mV peak, rise time less than 100 µs, 1ms pulse width, and 120 pulses / min repetition rate together with the ECG calibration signal CAL20000, and input to the ECG machine. The QRS curve of calibration signal can be clearly identified on ECG machine records, the difference between the recorded amplitude and the original amplitude of QRS curve does not exceed 20%, and the location of the pacemaker pulse can be clearly identified in the ECG machine records.

Appendix F EMC-Guidance and manufacture's declaration

A CAUTION

- The use of accessories, sensor and cable which exceed regulations can increase the electromagnetic emissions or decrease the electromagnetic immunity.
- Please don't use the equipment with tightness and pile. The equipment should be detected according to the requirements to make sure it working normally.
- Please make the equipment safety according to the EMC, and install, fix it under the situation required by EMC.
- Even though the other equipments which measure up the CISPR also may cause disturbance.
- When the amplitude of import signal narrow than regulated minimum amplitude may lead to inaccurate measure.
- Portable and mobile communications equipment will affect the performance of the equipment.

F.1 Guidance and manufacturer's declaration – electromagnetic emission

Guidance and manufacturer's declaration - electromagnetic emission

The Digital Electrocardiograph is intended for use in the electromagnetic environment specified below. The customer or the user of the ECG machine should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR11	Group 1	The ECG machine uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class A	
Harmonic emissions EN 61000-3-2	Class A	The ECG machine is suitable for use in all establishments, including domestic establishments and those directly connected to the public
Voltage fluctuations / flicker emissions EN 61000-3-3	Complies	low-voltage power supply network that supplies buildings used for domestic purposes.

F.2 Guidance and manufacture's declaration-electromagnetic immunity for all EQUIPMENT and SYSTEMS.

Guidance and manufacture's declaration-electromagnetic immunity

The ECG machine is intended for use in the electromagnetic environment specified below. The customer of the user of the ECG machine should assure that it is used in such an environment

Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment
			guidance
Electrostatic discharge	±6KV contact	±6KV contact	Floors should be wood, concrete or
(ESD)	±8KV air	±8KV air	ceramic tile. If floor are covered with
IEC61000-4-2			synthetic material, the relative
			humidity should be at least 30%
Electrical fast	±2KV for power supply lines	±1KV for power supply	Mains power quality should be that of
transient/burst	±1KV for input/output lines	lines	a typical commercial or hospital
IEC61000-4-4		±1KV for input/output lines	environment.
Surge IEC61000-4-5	±1KV differential mode	±1KV differential mode	Mains power quality should be that of
	±2KV common mode	±2KV common mode	a typical commercial or hospital
			environment.
Voltage dips, short	<5% UT	<5% UT	Main power quality should be that of
interruptions and	(>95% dip in UT)	for 0.5 cycle	a typical commercial or hospital
voltage variation power	for 0.5 cycle		environment. If the user of the ECG
supply input lines	40% UT	40% UT	machine requires
IEC61000-4-11	(60% dip in UT)	for 5 cycle	Continued operation during power
	for 5 cycle		mains interruptions, it is
	70% UT	70% UT	recommended that the ECG machine
	(30% dip in UT)	for 25 cycles	be powered from an uninterruptible
	for 25 cycles		power supply or a battery.
	<5% UT	<5% UT	
	(>95% dip in UT)	for 5 sec	
	for 5 sec		

^{-- 82 --} User Manual for Electrocardiograph

Power frequency	3A/m	3A/m	Mains power quality should be that of
(50Hz) magnetic			a typical commercial or hospital
field			environment.
IEC61000-4-8			
NOTE: UT is the A.C mains voltage prior to application of the test level.			

F.3 Guidance and manufacture's declaration-electromagnetic immunity for EQUIPMENT and SYSTEMS those are not LIFE-SUPPORTING

Guidance and manufacture's declaration-electromagnetic immunity

The ECG machine is intended for use in the electromagnetic environment specified below. The customer of the user of the ECG machine should assure that it is used in such an environment.

Immunity test	IEC60601 test level	Compliance	Electromagnetic environment-guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the ECG machine, including cables, than the recommended separation distance calculated from the
			equation applicable to the frequency of the transmitter. Recommended separation distance
Conducted RF IEC61000-4-6	3 Vrms 150kHz -80 MHz	3 Vrms	$d = \left[\frac{3.5}{V_1}\right] \sqrt{P}$
Radiated RF	3 V/m	3 V/m	$d = \left[\frac{3.5}{E_1}\right] \sqrt{P} \text{ 80MHz-800 MHz}$
IEC61000-4-3	80 MHz -2.5 GHz		$d = \left[\frac{7}{E_1}\right] \sqrt{P} \text{ 800MHz-2.5 GHz}$
			Where P is the maximum output power rating of the transmitter
			in watts (W) according to the transmitter manufacturer and d is
			the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, asdetermined by an

	electromagnetic site survey, range ^a should be less than the
	compliance level in each frequency range. ^b
	Interference may occur in the vicinity of Equipment marked with
	the following symbol:

NOTE1: At 80MHz and 800MHz, the higher frequency range applies.

NOTE2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ECG machine should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the ECG machine.
- b. Over the frequency range 150kHz to 80 MHz, field strengths should be less than 3V/m.

F.4 Recommended separation distance between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM for EQUIPMENT or SYSTEM that are not LIFE-SUPPORTING

Recommended separation distance between portable and mobile RF communications equipment and the ECG machine

The ECG machine is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ECG machine can help prevent electromagnetic Interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ECG machine as recommended below according to the maximum output power of the communications equipment.

	Separation distance according to frequency of transmitter (m)		
Rated maximum output	150 kHz-80 MHz	80 MHz-800 MHz	800 MHz-2.5 GHz
power of transmitter (W)	$d = \left[\frac{3.5}{V_1}\right] \sqrt{P}$	$d = \left[\frac{3.5}{E_1}\right] \sqrt{P}$	$d = \left[\frac{7}{E_1}\right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1: At 80MHz and 800MHz, the higher frequency range applies.

NOTE2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

-- Blank Page --

