

# Electronic Suction 1 Litre



**SS-1000**



Features	with a high flow, high vacuum and fashionable stylish		
	Adopting oil -free piston pump		
	Adjustable pressure button to select any safety pressure		
	Overflow-protection design		
	Bi-rich patented products		
Specifications	Power		AC220~230V 50Hz
	Input Power		160VA
	Reservoir volume		1000 ML/PC
	Noise level		≤60 dB(A) ( at 1 meter distance )
	Max negative pressure value		≥70Kpa ( 0.7 bar, 525mmHg )
	Range of negative pressure		10KPa~ max negative pressure
	Pumping Flow rate		≥16 L/MIN
	Operation environment		wozazhid
	Storage environment		Temperature: -20°C to 70°C ; Humidity: 15 to 93% RH
	Dimension		Length 280mm x Width 195 mm x Height 245mm
	Weight		3.0 KG
Accessories	Soft tube	Q'ty	3 pcs (2 pcs with 0.15M & 1 pcs with 1.3M)
	Suction tube	Q'ty	2 pcs (1 pcs for Adult & 1 pcs for Child)
	Bottle	Q'ty	1 pcs (1000ML)
	Filter	Q'ty	