

**Spirit<sup>®</sup>**  
SINCE 1975





# CK-A605T

Grandeur Series Lightweight Adult Dual Head Stethoscope

## Features:

- Lightweight chestpiece with an adult diaphragm for precision acoustics
- Binaural assembly with flexible outer spring
- Metal retaining ring



# CK-A604T

Identical features as CK-A605T but with pediatric chestpiece



Standard Color of PVC Y-Tubings



Standard Color of Soft Eartips



# CK-A603T

Grandeur Series Lightweight Nurse Single Head Stethoscope



# CK-608T

Majestic Series Neonatal Dual Head Stethoscope

## Features:

- Chrome plated with newborn type of chestpiece
- Binaural assembly with flexible outer spring
- Metal retaining ring



Standard Color of PVC Y-Tubings



Standard Color of Soft Eartips





# CK-S601PF

Deluxe Series Adult Dual Head Stethoscope

## Features:

- Stainless steel chestpiece provides excellent acoustic perception
- Binaural assembly with built-in sound centralized tubing
- Retaining and non-chill rings both in black or grey color
- Soft-sealing eartips in black or grey color
- Spare eartips and ID tag included



# CK-S606PF

Identical features as CK-S601PF but with pediatric chestpiece.



# CK-S607

Identical features as CK-S601PF but with infant chestpiece.



Diaphragm PF

Standard Color of PVC Y-Tubings



Standard Color of Soft Eartips



# CK-SS601PF

## Deluxe Series Dual Head Stethoscope

The revolution of interchangeable auscultation head that depends on diagnostic situations to use different convertible fittings.

### Features:

- CK-S601PF with one end fixed with adult diaphragm, the other end can be interchanged with an adult bell or child diaphragm
- Binaural assembly with built-in sound centralized tubing
- Retaining and non-chill rings both in black or grey color
- Soft-sealing eartips in black or grey color
- Spare eartips and ID tag included



# CK-625P

## Grandeur Series Teaching Stethoscope

### Features:

- 2 sets of binaural are applicable for 2 persons doing diagnostic proceedings simultaneously
- Binaural assembly with flexible outer spring



# CK-SS747PF

Deluxe Series Cardiology Stethoscope

## Features:

- Chestpieces and binaural are made of stainless steel
- Combination of a adult chestpiece and a convertible pediatric chestpiece for superior acoustic sensitivity
- Extra adult bell for low frequency response
- Binaural assembly with dual lumen tubing
- Retaining and non-chill rings both in black or grey color
- Soft-sealing eartips in black or grey color
- Spare eartips and ID tag included



# CK-S748PF

Deluxe Series Spirit III Single Head Cardiology Stethoscope

## Features:

- Stainless steel chestpiece with a tunable diaphragm for a superior acoustic sensitivity
- Binaural assembly with dual lumen tubing
- Retaining and non-chill rings both in black or grey color
- Soft-sealing eartips in black or grey color
- Spare eartips and ID tag included



Diaphragm PF

Standard Color of PVC Y-Tubings



01 02 03 05 07 10 16 21

Standard Color of Soft Eartips



01 02



# CK-110 CK-110L/F

Aneroid Portable Sphygmomanometer

## Specifications:

- Measurement Range: 0-300mmHg
- Inflation Method: Inflation and Air Release by Manual
- Cuff Size: approx. 510 x 142 mm / Nurse Type
- Cuff Material: Nylon

- \* Scale Glow-in-Dark
- \* Include a portable storage pouch.
- \* White reading scale is available



## CK-113

Palm Style Night Glow Aneroid Sphygmomanometer

## Specifications:

- Measurement Range: 0-300 mmHg
- Inflation Method: Inflation and Air Release by Manual
- Cuff Size: approx. 510 x 142 mm / Nurse Type
- Cuff Material: Deluxe Nylon or Cotton

- \* Include a portable storage pouch.





# CK-141L/F

## Wall Model Large Face Aneroid Sphygmomanometer

### Specifications:

- Measurement range: 0-300 mmHg
- Cuff size: approx. 510 x 142 mm / Nurse Type
- Cuff material: Nylon
- Adult size inflation system and coil tubing
- Accurate and reliable
- Easy to read large dial
- Scale Glow-in-Dark

\* Specialize: Wall hanged type, fits all sickbeds, dialysis room, operating room and other medical specialization room.



# CK-141AL/F

Features same as CK-141 but with cuff basket.



# CK-143

## Table Top Model

Identical features as CK-141 but with heavy base and save space.





# CK-142

Wall Model Large Face Aneroid Sphygmomanometer

## Mobile Type

Same as CK-141 but with

- Storage basket
- Wheel Color: Cream White or Black (CK-142C)
- Patented heavy wheel based, drag steady without falling
- Height adjustable: 870 mm (lowest), 1220 mm (highest)
- Wheel width: 410 mm



# CK-101

Mercurial Desk Model Sphygmomanometer

Specifications:

- Measurement Range: 0-300 mmHg
- Dimensions: 340 x 100 x 55mm
- Cuff Size: approx. 510 x 142 mm / Nurse Type
- Cuff Material: Nylon
- Product Color: Silver



# P-124

Vinyl Coated Wire Basket





# CK-907

Instrulite™

An innovative high lux and cool light diagnostic penlight, pocket clip press activator switch, pupil gauge and slide ruler. Comes complete with 2"AAA" batteries and individual package.



# CK-908

Instrulite™

## Radiantlite™

Made from a lightweight and sturdy metal case, features a secured high lux cool light bulb, "TipTop™" push-button switch and pocket clip. Comes in with (optional) or without imprinted pupil gauge/side ruler, and in different finishes, i.e. anodized, matte, or "Duracote™" solid colors.



# CK-909

Instrulite™

## Radiantlite™

This model offers a sturdy lightweight design, as well as a two-way on / off switch for convenient use. The "TipTop™" push-button switch allows the unit to be activated in short bursts. While the pocket clip can lock the unit in the "on" position. This alternative can extend battery life when usage.



*\*Customized paint coating is available, MOQ: 500 PCS/Color*



# CK-504

Queen Square Hammer

## Queen Square Hammer

A 4.65 cm Diameter chrome plated head with rubber ring and 30.5 cm long flexible handle with pointed tip.



## Rubber Bulb

### P-117

Standard Rubber Bulb with Valve.



### P-117/LF

Standard Latex Free Rubber Bulb with Valve.



### P-123

Long Coil Rubber Tubing

**Length:**

Normal: 31 cm

Stretched: 260 cm



### P-123L/F

Long Latex Free Coil Tubing

**Length:**

Normal: 31 cm

Stretched: 260 cm

**DISTRIBUTED BY**