

Spirit[®]
SINCE 1975





CK-A605T

Grandeur Series Lightweight Adult Dual Head Stethoscope

Features:

- Lightweight chestpiece with an adult diaphragm for precision acoustics
- Binaural assembly with flexible outer spring
- Metal retaining ring



CK-A604T

Identical features as CK-A605T but with pediatric chestpiece



CK-A603T

Grandeur Series Lightweight Nurse Single Head Stethoscope



Standard Color of PVC Y-Tubings



Standard Color of Soft Eartips



CK-608T

Majestic Series Neonatal Dual Head Stethoscope

Features:

- Chrome plated with newborn type of chestpiece
- Binaural assembly with flexible outer spring
- Metal retaining ring



Standard Color of PVC Y-Tubings



Standard Color of Soft Eartips





CK-S601PF

Deluxe Series Adult Dual Head Stethoscope

Features:

- Stainless steel chestpiece provides excellent acoustic perception
- Binaural assembly with built-in sound centralized tubing
- Retaining and non-chill rings both in black or grey color
- Soft-sealing eartips in black or grey color
- Spare eartips and ID tag included



CK-S606PF

Identical features as CK-S601PF but with pediatric chestpiece.



CK-S607P

Identical features as CK-S601PF but with infant chestpiece.



Diaphragm PF

Standard Color of PVC Y-Tubings



01 02 03 04 05 06 07 08 09 10 11 12 13 14 15 16 17 18



20 21

Standard Color of Soft Eartips



01 02

CK-SS601PF

Deluxe Series Dual Head Stethoscope

The revolution of interchangeable auscultation head that depends on diagnostic situations to use different convertible fittings.

Features:

- CK-S601PF with one end fixed with adult diaphragm, the other end can be interchanged with an adult bell or child diaphragm
- Binaural assembly with built-in sound centralized tubing
- Retaining and non-chill rings both in black or grey color
- Soft-sealing eartips in black or grey color
- Spare eartips and ID tag included



CK-625P

Grandeur Series Teaching Stethoscope

Features:

- 2 sets of binaural are applicable for 2 persons doing diagnostic proceedings simultaneously
- Binaural assembly with flexible outer spring



CK-SS747PF

Deluxe Series Cardiology Stethoscope

Features:

- Chestpieces and binaural are made of stainless steel
- Combination of a adult chestpiece and a convertible pediatric chestpiece for superior acoustic sensitivity
- Extra adult bell for low frequency response
- Binaural assembly with dual lumen tubing
- Retaining and non-chill rings both in black or grey color
- Soft-sealing eartips in black or grey color
- Spare eartips and ID tag included



CK-S748PF

Deluxe Series Spirit III Single Head Cardiology Stethoscope

Features:

- Stainless steel chestpiece with a tunable diaphragm for a superior acoustic sensitivity
- Binaural assembly with dual lumen tubing
- Retaining and non-chill rings both in black or grey color
- Soft-sealing eartips in black or grey color
- Spare eartips and ID tag included



 *Standard Color of PVC Y-Tubings* *Standard Color of Soft Eartips*

| | | | | | | | | | |
|---|---|---|---|---|---|---|---|---|---|
|  |  |  |  |  |  |  |  |  |  |
| 01 | 02 | 03 | 05 | 07 | 10 | 16 | 21 | 01 | 02 |

Diaphragm PF



CK-110

Aneroid Portable Sphygmomanometer

Specifications:

- Measurement Range: 0-300mmHg
- Inflation Method: Inflation and Air Release by Manual
- Cuff Size: approx. 510 x 142 mm / Nurse Type
- Cuff Material: Nylon

- * Scale Glow-in-Dark
- * Include a portable storage pouch.
- * White reading scale is available



CK-112P

Professional Switch Aneroid Sphygmomanometer

Specifications:

- Measurement Range: 0-300mmHg
- Product Rotate: 360°
- Inflation Method: Inflation and Air Release by Manual
- Cuff Size: approx. 510 x 142 mm / Nurse Type
- Cuff Material: Nylon

- * Include a portable storage pouch.
- * Gauge protector, gauge weight is optional.



CK-141

Wall Model Large Face Aneroid Sphygmomanometer

Specifications:

- Measurement range: 0-300 mmHg
- Cuff size: approx. 510 x 142 mm / Nurse Type
- Cuff material: Nylon
- Adult size inflation system and coil tubing
- Accurate and reliable
- Easy to read large dial
- Scale Glow-in-Dark

* Specialize: Wall hanged type, fits all sickbeds, dialysis room, operating room and other medical specialization room.



CK-141A

Features same as CK-141 but with cuff basket.



CK-143

Table Top Model

Identical features as CK-141 but with heavy base and save space.





CK-142

Wall Model Large Face Aneroid Sphygmomanometer

Mobile Type

Same as CK-141 but with

- Storage basket
- Wheel Color: Cream White or Black (CK-142C)
- Patented heavy wheel based, drag steady without falling
- Height adjustable: 870 mm (lowest), 1220 mm (highest)
- Wheel width: 410 mm



CK-101

Mercurial Desk Model Sphygmomanometer

Specifications:

- Measurement Range: 0-300 mmHg
- Dimensions: 340 x 100 x 55mm
- Cuff Size: approx. 510 x 142 mm / Nurse Type
- Cuff Material: Nylon
- Product Color: Silver



P-124

Vinyl Coated Wire Basket



CK-907

Radiantlite™

An innovative high lux and cool light diagnostic penlight, pocket clip press activator switch, pupil gauge and slide ruler. Comes complete with 2"AAA" batteries and individual package.



CK-908

Radiantlite™

Radiantlite™

Made from a lightweight and sturdy metal case, features a secured high lux cool light bulb, "TipTop™" push-button switch and pocket clip. Comes in with (optional) or without imprinted pupil gauge/side ruler, and in different finishes, i.e. anodized, matte, or "Duracote™" solid colors.



CK-909

Radiantlite™ II

Radiantlite™

This model offers a sturdy lightweight design, as well as a two-way on / off switch for convenient use. The "TipTop™" push-button switch allows the unit to be activated in short bursts. While the pocket clip can lock the unit in the "on" position. This alternative can extend battery life when usage.



**Customized paint coating is available, MOQ: 500 PCS/Color*



Hammer



CK-501

Taylor Percussion Hammer

Traditional style with 19.0cm. long chrome plated handle with triangular rubber head.



CK-504

Queen Square Hammer

A 4.65 cm Diameter chrome plated head with rubber ring and 30.5 cm long flexible handle with pointed tip.



CK-507

Babinski Telescoping Percussion Hammer

Extendable chrome plated hammer with pointed tip. Head (4.8 cm diameter) maybe used in both horizontal and vertical positions without unscrewing the handle.



ACCESSORIES

Eartip & Shield Name Tag Selection



P-212

Small Size Soft Rotating Eartip



P-213

Large Size Soft Rotating Eartip



P-251

Shield Name Tag

Cuff Selection



P-101NAD

Nylon Nurse / Adult Cuff with D-ring
52cm x 14.2cm



P-102NI

Nylon Infant Cuff
24.5cm x 8cm



P-102NA

Nylon Nurse / Adult Cuff
51cm x 14.2cm



P-102NP

Nylon Pediatric Cuff
35cm x 11cm



P-103 NDX

Deluxe Nylon Calibrated Cuff
51cm x 14.2cm



P-106 NXL

Nylon Obese / Large Cuff
67.5cm x 8cm

Bladder Selection



P-107 ADT

Adult 2-Tube Rubber Bladder
22.5cm x 13cm



P-107 ADT / LF

Latex-Free Adult 2-Tube Rubber Bladder
22.5cm x 13cm



P-107 AST

Adult Single Tube Rubber Bladder
22.5cm x 13cm



P-107 IDT

Infant 2-Tube Rubber Bladder
11.5cm x 7cm



P-107 IST

Infant Single Tube Rubber Bladder
11.5cm x 7cm



P-107PDT

Pediatric 2-Tube Rubber Bladder
17cm x 10cm



P-107PST

Pediatric Single Tube Rubber Bladder
17cm x 10cm



Rubber Bulb

P-117

Standard Rubber Bulb with Valve.



P-117/LF

Standard Latex Free Rubber Bulb with Valve.



P-123

Long Coil Rubber Tubing

Length:

Normal: 31 cm

Stretched: 260 cm



P-123L/F

Long Latex Free Coil Tubing

Length:

Normal: 31 cm

Stretched: 260 cm

DISTRIBUTED BY

