

**ΩMEGA**  
Health Care

# Oxygen Concentrator



**P-10**



- Altitude training allows athlete`s bodies to produce more red blood cells.
- Simulated Altitude Training is quickly emerging as the world's most promising development in the world of athletic performance and peak fitness.
- As your body begins to adjust to repeated exposure to reduced oxygen (hypoxic) environments, amazing physiological changes occur making you faster, stronger, healthier and leaner.

## Specification

Rated Flow	0-10L/min
Hypoxic Air Purity	9.5%-20.5%
Maximum Flow	120L/min
Power Consumption	≤ 530
Working Voltage	AC: 220/110 V ± 10% 50/60 Hz ±1
Altitude	0-8500m
Outlet Pressure	0.04-0.08Mpa
Noise Level	≤ 50dB
Net Weight	31kg
Gross Weight	33kg
Dimension	365X400X650mm(L*W*H)
Standard Function	Over Heat Alarm , Power Failure Alarm ,Timing Function , Working Hours Display.

**ΩMEGA**  
Health Care

# Oxygen generator with Nebulizer



**(VH5-N)**



#### Product Advanced Features:

- Waterproof structure of Oxygen Cylinder to keep longer lifespan.
- Power Modulation Technique to achieve stable working continuously.
- Self-check system provide light and sound alarm with data indicating on LCD when failure including low purity, power failure, overheat, over/low pressure, and control valve failure.
- Include Nebulizer.

#### Standard Functions

- Timing function for working hours control.
- 4 layer bacteria filter inside to keep bacteria away.
- Suitable for 24 hours continuous operation.
- LCD showing cumulative working hours and present working hours.

#### Condition For Transportation And Storage

- Environment temperature scale : -20°C to 45°C.
- Comparative humidity scale: ≤95%.
- Air pressure scale: 500hpa-1060hpa.

**Note: Please place the machine at a good ventilation area and make sure that the exhaust port underneath the machine is not covered.**

Model	(VH5-N)
Flow Rate	0-5 L/min
Power Consumption	350 W
Noise Level	≤40 dB
Purity	93±3%
Outlet Pressure	0.04-0.06 Mpa
Working Voltage	220V±10% 50Hz±1%;
Dimension	High 24.5" x wide 15.3" x deep 13.5"
Net Weight	21 kg
Gross Weight (carton)	24 kg
Warranty	1 year
Certificate	ISO13485 CE FDA
Working Hours	up to 10000 hours

Address : 3 Shenton Way #16-08 Shenton House Singapore 068805  
 Email: [international@omega-h-c.com](mailto:international@omega-h-c.com) [WWW.omega-h-c.com](http://WWW.omega-h-c.com)